BD Surgical Clipper and Charging Adapter Stand Instructions For Use

REF 5513E / 5514A / 5514E / 5514U / 5514K

en English 3

fr Français 15

de Deutsch 27

it Italiano 39

es Español 51

pt Português 63

nl Nederlands... 75

sv Svenska 87

da Dansk 99

fi Suomi 111


no Norsk 123






el Greece 135

tr Türkçe..... 147

pl Polski 159

cz Česky 171

Parts Identification (diagram #1):

1. REF 4406, 4403A or 4412A: Disposable blade, single use only
 REF 4406: General purpose surgical blade designed for all types of hair (light blue)
 REF 4403A: SensiClip™ surgical blade designed for clipping areas of loose or sensitive skin (light pink)
 REF 4412A: Neuro surgical blade designed for head hair and/or thick body hair (blue)
2. Clipper handle (REF 5513E)
3. On/Off switch
4. Recharge () LED
5. Charge status () LED
6. Charging adapter (REF 5514 series): CE marked except REF 5514A. (North America version) Be aware of appropriate plug configuration for your region.

Intended Use

The BD Surgical Clipper REF 5513E is a rechargeable clipper used with charging adapter (REF 5514 series) and BD disposable blades only (REF 4406, 4403A or 4412A). It is intended to remove head and body hair prior to any medical procedure requiring hair removal. The hair is removed by each blade oscillated by an electric motor. The BD 5513E Clipper will easily and effectively remove body hair and even the thickest hair from the chest with BD disposable blade REF 4406 BD, scalp and other thick coarse hair with BD disposable blade REF 4412A, and other difficult-to-clip areas of the body with BD disposable blade REF 4403A. The Clipper effectively removes wet or dry hair.

Charging

All Clippers arrive discharged. Before patient use, it is recommended to charge the Clipper a minimum of overnight (Maximum fully charge 17 hours). Switch Clipper off to charge. Place charger upright on a flat surface or wall mount. The charge status () LED glows to indicate it is receiving a charge. Except when charging, it is recommended to unplug the charging adapter from the electrical outlet. This will reduce overall energy consumption and is compliant with best practices. Below are the LED indicators on the Clipper:

-  Orange Solid. Battery is charging
-  Off Battery is fully charged
-  Blue Flashing. Battery charge low (5 minutes of runtime remaining)

Wall mount of charging adapter (diagram #2):

The charging adapter can be hung on the wall using wood screws. Place charging adapter per healthcare facility guidelines.

1. Select an appropriate location on the wall where the wall mount can be securely attached.
2. Attach the wall mount to the wall as shown with wood screws.
 - Use no. 6 by 1 inch (25.4mm) round head screws when fixing to wood; or no. 4-6 Molly or Anchor bolts when fixing to sheet rock.
3. Hook the charger onto these screws.
4. Plug the power cord into a power outlet.

Instructions for use

Only trained healthcare personnel should use the Surgical Clipper. If the Surgical Clipper is used for the first time, see the training materials provided by BD. Health care personnel should wear gloves when performing hair removal. Clipper should be disinfected prior to initial use.

Suitable for use in healthcare facility environment.

Attaching a disposable blade to the Clipper handle (diagram #3):

Make sure the Clipper is turned off. Using gloved hands, peel back and remove cover from the blade blister package. Position the new blade on top of the Clipper handle, aligning the mark on the side of the blade with the first line on the handle. Slide the blade back in the direction of the arrow until it clicks into place and the mark on the side of the blade aligns with the dot on the Clipper handle. If blade is dropped or becomes contaminated, discard it.

Patient's skin should be clean. Healthcare professional should instruct patient to avoid sudden movement during clipping process. Use dry or wet clipping method. Choose the appropriate blade depending on the area to clip. Attach blade (see above). Switch Clipper on (diagram #1 – mark 3). Check Clipper and blade for abnormal vibration or noise before patient use. See operating noise specifications listed in the Surgical Clipper specifications table. If abnormalities are found do not use the Clipper. With logo facing up, hold Clipper at approximate 40° angle to the patient, (hold between thumb and forefinger with hand underneath, similar to using a pen or a pencil, or with a standard over grip). Glide face of blade gently against surface of skin (diagram #4). For best results clip against the direction of hair growth. It is best to always stretch skin taut for clipping.

Removing blade:

After clipping, switch Clipper off to stop operation. Position the Clipper blade facing down (diagram #5) over an appropriate trash receptacle, and use your thumb to push the blade forward along the blade frame releasing and disposing of the used blade. Discard blade per facility protocol.

Disinfection:

Clipper should be disinfected after each use. The Clipper may be rinsed under running water or submerged, for up to 30 minutes, in water to remove any loose hair, debris or other gross contaminants. Cleaning the Clipper can be achieved with a pre-moistened wipe or cloth containing a sodium hypochlorite or hydrogen peroxide solution, following the manufacturer's instructions for use. Clipper should then be wiped with a pre-moistened disinfectant wipe or cloth moistened with a disinfectant solution. When necessary to clean the charging adapter, disconnect from wall outlet and wipe with the compatible disinfectants as specified for the Clipper. Compatible disinfectants include sodium hypochlorite solution, quaternary ammonium/isopropyl alcohol solution, hydrogen peroxide and benzalkonium chloride/dodecyl dimethyl ammonium chloride/PHMB solutions. Refer to disinfectant instructions for use. Allow Clipper to dry completely before returning to the charging adapter. Disinfectants containing benzyl alcohol can cause chemical damage to the external housing and should not be used for routine cleaning. Use of benzyl alcohol will void the warranty period stated for this product. Should the Clipper become grossly contaminated, discard and contact your local sales representative for a replacement.

WARNING:

BD blades (REF 4406, 4403A, and 4412A) are single use only and specifically designed for use with the BD 5513E. The user assumes responsibility for appropriate use of this Clipper. Using blades not manufactured or approved by BD will void any warranty and patient results cannot be predicted. Re-use of blades may result in a nonfunctional product and could contribute to cross contamination; potentially putting patient safety at risk.

Surgical Clipper Specifications:

Mode of Operation	Continuous
BD Disposable Blades Mode of Operation	Single Use
Operating Voltage / Current	3.25V DC / 390mA (max.)
Weight	121 grams (4.3 oz.)
Dimension	6.0" (153 mm) H x 1.9" (49 mm) W x 1.5" (39 mm) D
Full Charge Time	Maximum 17 hours
Continuous Run Time (Fully Charged)	Disposable Blades REF 4406 & REF 4403A: Minimum of 80 minutes continuous run-time up to 120 minutes. Disposable Blade REF 4412A: Minimum of 50 minutes continuous run-time.
Operating Noise	REF 4406 Blade is less than 72dB. REF 4412A Blade is less than 69dB. REF 4403A Blade is less than 75dB.
Type of Protection Against Electric Shock	Internally Powered
Degree of Protection Against Electric Shock	Type B Applied Part
Degree of Protection Against Flammable Anesthetic	Equipment not suitable for use in the presence of a flammable anesthetic.
Ambient Operating Temperature, Humidity, and Atmospheric Pressure	+41 °F (+5 °C) to +95 °F (+35 °C), 30% to 75% RH, 70 to 106 kPa.
Storage Temperature, Humidity, and Atmospheric Pressure	+14 °F (-10 °C) to +140 °F (+60 °C), 10% to 100% RH, 50 to 106 kPa.

Surgical Clipper Specifications (continued):

Waterproofing	IPX7 Rating
Battery Type	Li-ion 680mAh
Battery Life	Minimum 350 Cycles

Charging Adapter Specifications:

Weight	200-250g ± 10g. (7-8.80 Oz.)
Power Source	100-240 V AC; 50-60 Hz, 30mA
Type of Protection Against Electric Shock	Class II
Dimension	3.15" (80 mm) H x 2.40" (61 mm) W x 1.73" (44 mm) D & Cord Length 73.62" (1870 mm)
Ambient Operating Temperature, Humidity, and Atmospheric Pressure	+41 °F (+5 °C) to +95 °F (+35 °C), 30% to 75% RH, 70 to 106 kPa.
Storage Temperature, Humidity, and Atmospheric Pressure	+14 °F (-10 °C) to +140 °F (+60 °C), 10% to 100% RH, 50 to 106 kPa.
Waterproofing	IPX7 Rating
Cord Plug Type	5514A - Type B 5514E - Type F 5514K - Type I 5514U - Type G

Symbol Explanation:

 Catalog Number	 Direct current	 Authorized representative in the European Community	 106 kPa 50 kPa Atmospheric pressure limitation
 Batch Code	 Watertight Equipment	 Class II (Double Insulation)	 100% 10% Humidity limitation
 Medical Device	 Alternating current	 Follow operating instructions	 +60°C -10°C Temperature limit
 Unique Device Identifier	 On / Off	 Do not reuse	 Type B applied part
 Manufacturer	 Use By	 Operating instructions	

Information for Users on Collection and Disposal of Old Equipment and used Batteries

These symbols on the products, packaging, and/or accompanying documents mean that used electrical and electronic products and batteries should not be mixed with general household waste. For proper treatment, recovery and recycling of old products and used batteries, please take them to applicable collection points, in accordance with your national legislation and the Directives 2002/96/EC and 2006/66/EC. By disposing of these products and batteries correctly, you will help to save valuable resources and prevent any potential negative effects on human health and the environment which could otherwise arise from inappropriate waste handling.

For more information about collection and recycling of old products and batteries, please contact your local municipality, your waste disposal service or the point of sale where you purchased the items. Penalties may be applicable for incorrect disposal of this waste, in accordance with national legislation.

Information on Disposal in other Countries outside the European Union

These symbols are only valid in the European Union. If you wish to discard these items, please contact your local authorities or dealer and ask for the correct method of disposal.

Caution

- The clipper REF 5513E must be used with REF 5514 series charger.
- Disposable single use blades (REF 4406, 4403A and 4412A) are designed for optimal use with 5513E only. Using blades not manufactured or approved by BD will void any warranty and patient results cannot be predicted.
- Do not submerge Clipper in water or other solution deeper than 3.3 feet (1 m) for no longer than 30 minutes.
- Do not plug in with wet hands. Do not connect to charging adapter if clipper is wet.
- Keep charging adapter cord away from heated surfaces.
- Do not position charging adapter where it is difficult to unplug.
- Do not place the Clipper on the charging adapter until charging adapter is seated on a flat surface or securely mounted to a wall.
- Prior to charging Clipper ensure charging adapter is free of metallic debris.
- Do not use near flammable anesthetic, aerosol spray or oxygen — administering equipment other than nasal or mask types.
- Do not expose to hot water, salt water, organic solvents or bleach solutions.
- Do not use hydrocarbon or phenol based cleaners or cleaners containing acetone or ketones.
- Do not use with damaged blade, handle, or both.
- Do not replace supply cord. If supply cord is damaged return to manufacturer or service agent for replacement, to avoid hazard.
- Do not sterilize.
- Do not take the housing apart as this can affect the watertight construction.
- Operating Clipper without blade could lead to injury.
- Service or repair must be performed only by qualified repair personnel. Service or repair performed by unqualified personnel may result in a risk of injury, electric shock or fire.

Caution (continued)

- Except when charging, it is recommended to unplug the charging adapter from the electrical outlet. This will reduce overall energy consumption and is compliant with best practices.
- Do not leave Clipper running without applying to skin for more than 1 minute as blade temperature may exceed 60 °C and potentially leading to thermal injury.
- During use, do not keep the Clipper blade applied to the same position of the patient's skin for longer than 1 minute (these operations may result in blade surface becoming hot).
- Inspect treatment site for selection of appropriate blade. If skin irregularities are present, proceed with caution.
- In cases of minor injury, seek medical treatment if necessary.
- Healthcare personnel should instruct patient to avoid sudden movement during clipping process.

Notes

- Avoid charging in direct sunlight or near radiating heat.
- When not used for extended periods of time, charge every six months to preserve battery life.
- Clipper contains a Lithium-ion battery. Dispose of at an officially assigned location, if available.
- The battery in this Clipper is not intended to be replaced by users. Diagram #6 is intended only for removal of the chargeable battery for the purpose of proper disposal.
- Do not apply Clipper blade to injured skin area.

Any serious incident that has occurred in relation to the device should be reported to BD at bd.com and the Competent Authority of the Member State.

For environmental protection and recycling of materials

This Clipper contains a lithium-ion battery. Remove the battery before disposal of the Clipper (#6). When discarding the battery and/or electrical equipment, it should not be mixed with general waste; it must be recycled or disposed of properly according to local regulations.

ATTENTION:

For the United States: A lithium ion battery that is recyclable powers the product you have purchased. Please call 1-800-8-BATTERY for information on how to recycle this battery.


Limited Warranty:

This Surgical Clipper product from BD, REF 5513E and REF 5514 series, is warranted against defects in material and workmanship under normal use for a period of 2 years from the date of purchase. If this product is determined to be defective within the warranty period, it will be repaired or replaced without charge by BD if you contact your local representative or customer service representative who will provide you with the proper authorization, necessary shipping labels and packing list. Please include a description of the defect with the returned unit. This warranty does not cover the replaceable blade unit nor does it cover damage caused by misuse, negligence, accident, or abuse, alteration or modification of the units by any other than BD or its authorized agent, or affixation of any accessory not manufactured or authorized by BD, or use of the Surgical Clipper with current or voltage other than that specified in these instructions.

The purchaser assumes responsibility for the safe operation of this Clipper in conjunction with approved and manufactured BD blades (REF 4406, 4403A, and 4412A). Only trained health care professionals or other appropriately trained individuals should use this equipment.

BD is not liable in either tort or contract for any cost, expense, loss, or damage whether direct, consequential, special, collateral, or incidental arising out of the use, misuse, or inability to use this product. Except as expressly provided above, BD disclaims all warranties including, but not limited to the implied warranties of merchantability and fitness for particular purpose.

#6

This equipment has been tested and found to comply with the limits for the medical devices to the IEC 60601-1-2, Medical Device Regulation 2017/745/EU. These limits are designed to provide reasonable protection against harmful interference. This equipment generates, uses and can radiate radio frequency energy and, if not installed and used in accordance with the instructions, may cause harmful interference to other devices in the vicinity. However, there is no guarantee that interference will not occur in a particular installation. If this equipment does cause harmful interference to the other devices which can be determined by turning the equipment off and on, the user is encouraged to try to correct the interference by one or more of the following measures;

- Reorient or relocate the receiving device.
- Increase the separation between the equipment.
- Connect the equipment into an outlet on a circuit different from that to which the other device(s) are connected.
- Consult the manufacturer or field service for help.

EMC Information of REF 5513E/5514 series

REF 5513E/5514 series needs special precautions regarding EMC and needs to be used according to the following EMC information.

WARNING:

- Use of REF 5513E/5514 series adjacent to or stacked with other equipment should be avoided because it could result in improper operation. If such use is necessary, REF 5513E/5514 series and the other equipment should be observed to verify that they are operating normally.
- Use of accessories, other than those specified or provided by the manufacturer of REF 5513E/5514 series could result in increased electromagnetic emissions or decreased electromagnetic immunity of REF 5513E/5514 series and result in improper operation.
- Portable RF communications equipment (including peripherals such as antenna cables and external antennas) should be used no closer than 30 cm (12 inches) to any part of REF 5513E/5514 series, including cables specified by the manufacturer. Otherwise, degradation of the performance of REF 5513E/5514 series could result.

Guidance and manufacturer's declaration — electromagnetic emissions

REF 5513E/5514 series is intended for use in the electromagnetic environment specified below.
The customer or the user of REF 5513E/5514 series should insure that it is used in such an environment.

Emission test	Compliance	Electromagnetic environment - guidance
RF emissions CISPR 11	Group 1	REF 5513E/5514 series uses RF energy only for its internal function. Therefore, its RF emissions are very low and are not likely to cause any interference in nearby electronic equipment. REF 5513E/5514 series is suitable use in all establishments, including domestic establishments and those directly connected to the public low-voltage power supply network that supplies buildings used for domestic purposes.
RF emissions CISPR 11	Class B	
Harmonic emissions IEC 61000-3-2	Class A	
Voltage fluctuations/ flicker emissions IEC 61000-3-3	Complies	

Guidance and manufacturer's declaration — electromagnetic immunity

REF 5513E/5514 series is intended for use in the electromagnetic environment specified below.
The customer or the user of REF 5513E/5514 series should insure that it is used in such an environment.

Immunity test	IEC 60601 test level	Compliance level	Electromagnetic environment - guidance
Electrostatic discharge (ESD) IEC 61000-4-2	±8 kV contact ±2/4/8/15 kV air	±8 kV contact ±2/4/8/15 kV air	Floors should be wood, concrete or ceramic tile. If floors are covered with synthetic material, the relative humidity should be at least 30%.
Electrical fast transient/burst IEC 61000-4-4	±2 kV for power supply lines 100kHz repetition frequency	±2 kV for power supply lines 100kHz repetition frequency	Mains power quality should be that of a typical commercial or hospital environment.
Surge IEC 61000-4-5	±0.5/1kV line to line	±0.5/1kV line to line	Mains power quality should be that of a typical commercial or hospital environment.
Voltage dips, short interruptions and voltage variations on power supply input lines IEC 61000-4-11	0% U_r ; 0.5 cycle at 0°, 45°, 90°, 135°, 180°, 225°, 270° and 315° 0% U_r ; 1 cycle and 70% U_r ; 25/30 cycles Single phase: at 0° 0% U_r ; 250/300 cycle	0% U_r ; 0.5 cycle at 0°, 45°, 90°, 135°, 180°, 225°, 270° and 315° 0% U_r ; 1 cycle and 70% U_r ; 25 cycles Single phase: at 0° 0% U_r ; 250 cycle	Mains power quality should be that of a typical commercial or hospital environment.
Power frequency (50/60 Hz) magnetic field IEC 61000-4-8	30 A/m 50 Hz/60Hz	30 A/m 50 Hz/60Hz	Power frequency magnetic fields should be at levels characteristic of a typical location in a typical commercial or hospital environment.

Note: U_r is the a.c. mains voltage prior to application of the test level

Guidance and manufacturer's declaration — electromagnetic immunity

REF 5513E/5514 series is intended for use in the electromagnetic environment specified below.

The customer or the user of REF 5513E/5514 series should insure that it is used in such an environment.

Immunity test	IEC 60601 test level	Compliance level	Electromagnetic environment - guidance
Conducted RF IEC 61000-4-6	3V 0.15MHz-80MHz 6V in ISM bands between 0.15MHz and 80MHz 80% AM at 1 kHz	3V 0.15MHz-80MHz 6V in ISM and amateur radio bands between 0.15MHz and 80MHz 80% AM at 1 kHz	Portable and mobile RF communications equipment should be used no closer to any part of REF 5513E/5514 series, including cables, than the recommended separation distance calculated from the equation applicable to the frequency of the transmitter. Recommended separation distance $d = 1.2 \sqrt{P}$ $d = 0.35\sqrt{P}$ 80MHz to 800MHz $d = 0.7\sqrt{P}$ 800MHz to 2.7GHz where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer and d is the recommended separation distance in meters (m). Field strength from fixed RF transmitters, as determined by an electromagnetic site survey, ^a should be less than the compliance level in each frequency range. ^a Interference may occur in the vicinity of equipment marked with the following symbol: 
Radiated RF IEC 61000-4-3	3V/m 80MHz-2,7GHz 80% AM at 1 kHz	10V/m 80MHz-2,7GHz 80% AM at 1 kHz	
Proximity fields from RF wireless communications equipment IEC 61000-4-3	IEC 60601-2: 2014 Table 9	IEC 60601-1-2: 2014 Table 9	$d = 6/\sqrt{fP}$ 380MHz to 5.8GHz

NOTE 1: At 80 MHz and 800 MHz, the higher frequency range applies.

NOTE 2: These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.

- a** Field strengths from fixed transmitters, such as base stations for radio (cellular/cordless) telephones and land mobile radios, amateur radio, AM and FM radio broadcast and TV broadcast cannot be predicted theoretically with accuracy. To assess the electromagnetic environment due to fixed RF transmitters, an electromagnetic site survey should be considered. If the measured field strength in the location in which REF 5513E/5514 series is used exceeds the applicable RF compliance level above, REF 5513E/5514 series should be observed to verify normal operation. If abnormal performance is observed, additional measures may be necessary, such as reorienting or relocating REF 5513E/5514 series.
- b** Over the frequency range 150 kHz to 80 MHz, field strengths should be less than 3 V/m.

Recommended separation distance between portable and mobile RF communications equipment and REF 5513E/5514 series

REF 5513E/5514 series is intended for use in an electromagnetic environment in which radiated RF disturbances are controlled. The customer or the user of REF 5513E/5514 series can help prevent electromagnetic interference by maintaining a minimum distance between portable and mobile RF communications equipment (transmitters) and REF 5513E/5514 series as recommended below, according to the maximum output power of the communications equipment.

Rated maximum output power of transmitter W	Separation distance according to frequency of transmitter m		
	150 kHz to 80 MHz $d = 1.2 \sqrt{P}$	80 MHz to 800 MHz $d = 0.35 \sqrt{P}$	800 MHz to 2.7 GHz $d = 0.7 \sqrt{P}$
0.01	0.12	0.035	0.07
0.1	0.38	0.11	0.22
1	1.2	0.35	0.7
10	3.8	1.1	2.2
100	12	3.5	7

For transmitters rated at a maximum output power not listed above, the recommended separation distance d in meters (m) can be estimated using the equation applicable to the frequency of the transmitter, where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer.

NOTE 1: At 80 MHz and 800 MHz, the separation distance for the higher frequency range applies.

NOTE 2: These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.

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Made in Japan
36-10431 Rev. F

Becton, Dickinson and Company
75 North Fairway Drive
Vernon Hills, IL 60061
USA

BD Switzerland Sàrl
Terre Bonne Park – A4
Route de Crassier 17
1262 Eysins, Switzerland

BD, the BD Logo, and SensiClip are trademarks of
Becton, Dickinson and Company or its affiliates.
©2020 BD. All rights reserved.

Becton Dickinson Pty Ltd.
4 Research Park Drive
North Ryde, NSW, 2113, Australia

Becton Dickinson Limited
14B George Bourke Drive
Mt. Wellington, Auckland, 1060 New Zealand

Sponsored in Australia and New Zealand by BD.

bd.com
bd.com/surgeryifu
bd.com/symbols-glossary

Issuance Date: 2020.02.05