Valleylab™

Smoke Evacuation Pencil

- en** Instructions for Use
- fr** Mode d'emploi
- pt** Instruções de uso
- da** Brugervejledning
- de** Bedienungsanleitung
- es** Instrucciones de uso
- fi** Käyttöohjeet
- it** Istruzioni per l'uso
- nl** Gebruksinstructies
- no** Bruksanvisning
- sv** Bruksanvisning
- ru** Инструкции по эксплуатации
- zh** 使用说明

Part No. PT00134829

STERILE EO

Caution, consult accompanying documents

Single use

COVIDIEN, COVIDIEN with logo, and Covidien logo and Positive Results for Life are U.S. and internationally registered trademarks of Covidien AG. Other brands are trademarks of a Covidien company,™ brands are trademarks of their respective owner.

May be covered by U.S. Patents:
www.covidien.com/patents

© 2013 Covidien.
Manufactured for Covidien Inc,
15 Hampshire Street, Mansfield, MA 02048 USA.
EC REP Covidien Ireland Limited,
IDA Business & Technology Park, Tullamore,
Ireland.

Made in China. Printed in China.

www.covidien.com

Consult instruction for use

REV 04/2021

REF CVPLP2000 **Smoke Evacuation Pencil**

This product contains DEHP. The intended use limits exposure to transient contact, minimizing the risk of DEHP release from the device. In order to avoid undue risk of DEHP exposure in children and nursing or pregnant women, product should only be used as directed.

Do not use if package is opened or damaged.

Not made with natural rubber latex.

For use with maximum voltage of 10 kV p-p.

This pencil is compatible with Covidien ES electrodes having a 3/32 inch (2.38 mm) diameter.

Operating Parameters: 10 to 40 °C (50 to 104 °F) at 30 to 75% relative humidity

Storage and Transport Parameters: -29 to 60 °C (-20 to 140 °F) at 30 to 85% relative humidity

Indications for Use

The Smoke Evacuation Pencil is designed for general electrosurgical applications, including cutting and coagulation, and for removing smoke generated by electrosurgery when used in conjunction with an effective smoke evacuation system. The pencil enables the operator to remotely conduct an electrosurgical current from the output connector of an electrosurgical unit to the target tissue for the desired surgical effect.

The Indications for use for the integrated Smoke Evacuation Pencil Include:

- To remove smoke plume from the surgical site.
- To remotely conduct an electrosurgical current from the output connector of an electrosurgical unit to the target tissue for the desired surgical effect.

Precaution

Sterility is guaranteed unless packaging is damaged or any seal is broken.

Precaution

The smoke evacuation pencil plug includes a standard three-prong plug that fits into current electrosurgical generators. The device is intended for use with electrosurgical generators. Refer to generator documentation to ensure compatibility.

The Valleylab™ Smoke Evacuation Pencil is not a fluid-removal device, therefore, it should not be used for such applications.

Warning

Do not bend or modify the electrode tip. Tip modifications may result in breakage or other damage. If the electrode is damaged, discard it.

Not all electrodes are designed for use in coagulation mode.

The CVPLP2000 contains di-(2-ethylhexyl) phthalate (DEHP) in concentrations above 0.1% by weight.

If the original electrode is removed from the pencil, visually confirm the new electrode is fully inserted and secured before activating the pencil.

Do not remove the active electrode from the surgical site while energy is activated. Unintended tissue effect or burns may occur.

For procedures where visualization may be impaired, be alert to these potential hazards:

- The electrode tip may remain hot enough to cause burns after the current has been deactivated.
- Inadvertent activation or movement of the activated electrode outside of the field of vision may result in injury to the patient.
- Localized burns to the patient or physician may result from electrical currents carried through conductive objects. Current may be generated in conductive objects by direct contact with the active electrode or by the active accessory being in close proximity to the conductive object.

Do not use electrosurgery in the presence of flammable anesthetics or other flammable gases, fluids, or objects, or in the presence of oxidizing agents, as fire may result.

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
Warning

This device is for use only by trained, licensed physicians. Do not use electrosurgical equipment unless properly trained to use it in the specific procedure being undertaken. Use of this device without such training can result in serious, unintended patient injury. While every attempt has been made to reduce patient and user risks, all surgeries using this device carry some residual risk, even when used by trained physicians.

Do not use in patients who have electronic implants such as cardiac pacemakers without first consulting a qualified professional (for example, cardiologist). A possible hazard exists because interference with the action of the electronic implant may occur, or the implant may be damaged.

Confirm proper power settings prior to and during a procedure. Use the lowest power settings that achieve the desired effect. If increased power settings are requested, check the patient return electrode and all accessory connections before major power-setting adjustments.

Valleylab™ Smoke Evacuation Pencil

- ① Blade
- ② Smoke evacuation tube
- ③ Cut button
- ④ Coag button
- ⑤ 3-prong connector

Warning

Do not use with incompatible accessories. This pencil is designed to fit with Covidien electrodes having a 3/32 inch (2.38 mm) diameter.

Facial and other body hair is flammable. Water soluble surgical lubricating jelly may be used to cover hair close to the surgical site to decrease flammability.

Warning

Do not modify or add to the insulation of active electrodes. If insulated electrodes are required, use an appropriate insulated electrode.

Do not activate the pencil while installing or removing the electrode. The user may also turn the generator off, place the generator in lowest power settings, or unplug the pencil from the generator during manipulation of the electrode.

The electrode tip is hot. The surface of the active electrode may remain hot enough to cause burns after the RF current is deactivated.

Electric Shock Hazard: Do not connect wet accessories to the generator.

Do not contact or arc the electrode to metal surfaces.

Fire/Explosion Hazard: Sparking and heating associated with electrosurgery can provide an ignition source.

Fire/Explosion Hazard: Position surgical electrode cables to avoid contact with the patient or other leads.

Confirm proper electrosurgical generator setting before proceeding with surgery. Use the lowest power setting to achieve the desired effect.

Explosion Hazard: Do not use electrosurgery in the presence of flammable anesthetics.

Fire Hazard: Place active accessories in a holster or in a clean, dry, non-conductive and highly visible area away from the patient when not in use. Inadvertent contact with the patient may result in burns. Contact with drapes or linens may cause a fire.

Always use the lowest power setting that achieves the desired surgical effect. Use the active electrode for the minimum time necessary in order to reduce the possibility of unintended burn injury.

Using coated electrodes at high power settings may cause damage to the coating. If the coating is damaged, discard the electrode. The electrode coating may deteriorate when used with tissue response generators at higher power settings.

Warning

Fire Hazard: During Oropharyngeal Surgery Verify endotracheal tubes are leak free and that the cuff seals properly to prevent oxygen leaks.

If an uncuffed tube is in use, pack the throat with wet sponges around the uncuffed tube, and be sure to keep sponges wet throughout the procedure.

Question the need for 100% oxygen during oropharyngeal or head and neck surgery.

If necessary, scavenge excess oxygen with separate suction.

If possible, stop supplemental oxygen at least one minute before and during use of electrosurgery.

Some surgeons may elect to “buzz the hemostat” during surgical procedures. It is not recommended, and the hazards of such a practice probably cannot be eliminated. Burns to the surgeon’s hands may result. To minimize the risk, take these precautions:

- Do not “buzz the hemostat” with a needle electrode.
- Do not lean on the patient, the table, or the retractors while buzzing the hemostat.
- Activate CUT rather than COAG. CUT has a lower voltage than COAG.
- Use the lowest power setting possible for the minimum time necessary to achieve hemostasis.
- Activate the generator after the accessory makes contact with the hemostat. Do not arc to the hemostat.
- Firmly grasp as much of the hemostat as possible before activating the generator. This disperses the current over a larger area and minimizes the current concentration at the fingertips.
- “Buzz the hemostat” below hand level (as close as possible to the patient) to reduce the opportunity for current to follow alternate paths through the surgeon’s hands.
- When using a coated or nonstick blade electrode, place the edge of the electrode against the hemostat or other metal instrument.

Inspect accessories and cords for breaks, cracks, nicks, or other damage before every use. If damaged, do not use. Failure to observe this caution may result in injury or electrical shock to the patient or surgical team.

Visual inspection alone may not be sufficient to ensure that the insulation is intact.

Verify all anesthesia circuit connections are leak free before and during use of electrosurgery.

Warning

Electrosurgical accessories that are activated or hot from use can cause unintended burns to the patient or surgical personnel:

- Electrosurgical accessories may cause fire or burns if placed close to or in contact with flammable materials, such as gauze or surgical drapes. Place longer electrodes, such as extended electrodes, away from the patient and drapes.

Tissue buildup (eschar) on the tip of an active electrode may create embers that pose a fire hazard, especially in oxygen-enriched environments. Keep the electrode clean and free of all debris.

Conductive fluids (e.g., blood or saline) in direct contact with an active electrode or in close proximity to any active accessory may carry electrical current and cause unintended burns to the patient. This can happen as a result of either direct coupling with the active electrode or capacitive coupling between the active electrode and the external surface of the electrode insulation. Therefore, to prevent unintended burns in the presence of conductive fluids.

- Always keep the external surface of the active electrode away from adjacent tissue while activating the generator.
- Clear conductive fluid from the electrode before activating the electrosurgical device.
- Keep gauze and sponges wet.
- Keep electrosurgical electrodes away from flammable materials and environments enriched with oxygen (O_2).
- Use of electrosurgery in oxygen-enriched environments increases the risk of fire. Therefore, take measures to reduce the oxygen concentration at the surgical site. If possible, stop supplemental oxygen at least one minute before and during use of electrosurgery.
- Avoid enriched oxygen (O_2) and nitrous oxide (N_2O) atmospheres near the surgical site, especially during head and neck surgery. Both oxygen and nitrous oxide support combustion and may result in fires and burns to patients or surgical personnel.

- Prevent pooling of flammable fluids and the accumulation of flammable or oxidizing gases or vapors under surgical drapes or near the surgical site.
- Do not activate the generator until flammable vapors from skin prep solutions and tinctures have dissipated.
- Avoid the accumulation of naturally occurring flammable gases in body cavities such as the bowel.

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Warning

Due to concerns about the carcinogenic and infectious potential of electrosurgical by products (such as tissue smoke plume and aerosols), protective eye wear, filtration masks, and effective smoke evacuation equipment should be used in both open and laparoscopic procedures.

Pediatric applications and/or procedures performed on small anatomic structures may require reduced power settings. The higher the current flow and the longer the current is applied, the greater the possibility of unintended thermal damage to tissue, especially during use on small appendages. Consult the hospital biomedical engineer for further information.

Warning

Confirm that the active electrode is connected to the active generator port and for monopolar use ensure that the patient electrode is appropriate, properly applied to the patient, and connected to the generator.

To reduce the risk of an inadvertent electrosurgical burn at the monitoring electrode or probe site, place the electrode and/or probe as far away as possible from the electrosurgical site and/or patient return electrode. Protective impedances (resistors or RF inductors) installed in the monitoring leads may reduce the risk of such burns. Consult the hospital biomedical engineer for further information.


1. Open product packaging and remove the Valleylab™ Smoke Evacuation Pencil assembly. Be careful not to allow the device or any of its attachments to come in contact with the non-sterile portions of the packaging.

2. Connect the tubing provided with the Smoke Evacuation Pencil to the smoke evacuation system or other appropriate vacuum source.

3. **Prior to connecting the Smoke Evaporation Pencil to the electrosurgical generator,** press gently on the plastic shield covering the electrode to ensure the electrode is fully inserted.

4. Remove and discard the electrode's plastic shield.

5. Plug the pencil's 3-prong power cord connector into an electrosurgical generator's monopolar accessory port. Ensure the port's receptacles match the plug's pin configuration.

Electrosurgical Generators and Pencils

- Check the manufacturer's instructions for proper setup, use, and troubleshooting of the electrosurgical generator. Refer to the manufacturer's precautions before use.
- An adapter may be needed to connect the pencil to generators not manufactured by Covidien. Contact your Covidien sales representative for assistance in selecting the appropriate adapter.

Smoke Evacuation

- Check the manufacturer's instructions for proper setup, use, and troubleshooting of the smoke evacuation system. Refer to the manufacturer's precautions before use.

Before Surgery

Warning

Read all warnings, cautions, and instructions provided with the generator and other accessories used with this device



Please pay attention to further warnings in the operating instructions for the generator.

Before use, examine the generator and accessories for defects. Do not use cables or accessories with damaged insulation or connectors. If the accessories are damaged, discard them.

The electrode must fit completely and securely into the device. An incorrectly seated electrode may result in burns to the patient or surgical personnel.


Hold the pencil safely away from the patients, personnel, and surgical drapes during testing.

6. Attach the holster to the surgical drape by putting the drape material through the slots.

7. With the electrosurgical generator off, maximize smoke evacuation by positioning the smoke evacuation tube close to the point of the tissue interaction. Do not obstruct the view of the target tissue.

8. Follow the instructions for applying the return electrode to the patient and connecting the return electrode to the generator.

9. Test the pencil before surgery:

- Turn on the generator.
- Set the generator to 1 watt for both cut and coagulate modes.
- Press the yellow cut button on the pencil. Verify that the yellow cut indicator on the generator illuminates.
- Press the blue coag button on the pencil. Verify that the blue coag indicator on the generator illuminates.

10. Confirm that all power settings on the generator are appropriate for the procedure to be performed. Verify that the generator output does not exceed rated accessory voltage (both return electrode and pencil electrode).

Precaution

Always confirm correct power settings if the electrode in this device is replaced or exchanged.

The available high-power settings may cause patient injury or product damage. Verify that generator power settings are appropriate for the procedure, device, and return electrode being used.

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Do not use needles as monitoring electrodes during electrosurgical procedures. Inadvertent electrosurgical burns may result.

Precaution

Pediatric applications and/or procedures performed on small anatomic structures may require reduced power settings. The higher the current flow and the longer the current is applied, the greater the possibility of unintended thermal damage to tissue, especially during use on small appendages.

Needle electrodes are fragile. Handle them with care to avoid damage to the needle and injury to hospital personnel.

Activate the electrosurgical unit only when you are ready to deliver electrosurgical current and the active tip is in view near target tissue (especially if looking through an endoscope).

Deactivate the electrosurgical unit before the tip leaves the surgical site.

Some electrodes are not designed for use in the coagulation mode.

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Precaution

Place active accessories in a holster or in a clean, dry, non-conductive and highly visible area away from the patient when not in use. Inadvertent contact with the patient may result in burns. Contact with drapes or linens may cause a fire.

Do not exceed maximum power limits as stated in these or other instructions for use. Exceeding recommended power settings may result in patient injury or product damage.

To prevent unintended injury to patient or operating room staff personnel, the conductive metal post of the electrode must be properly and completely installed into the nose of the pencil. Any exposed metal can result in arcing at the electrode and pencil connection.

Do not activate electrodes while in contact with or near other instruments, including cannulas. Localized burns to the patient or physician may result.

Do not activate the generator in an open circuit condition. To reduce the chances of unintended burns, activate the generator only when the active electrode is near or touching the target tissue.

Activate the generator only when ready to deliver electrosurgical current and when the active tip is in view, and near target tissue.

Notice

Needle electrodes are designed for precise low power use during monopolar electrosurgery. Using a needle at high power settings for extended periods of time may result in damage to the needle. Use low power settings for short periods of time to prevent damaging the needle tip.

When using an electrode, inspect the electrode frequently to verify that the coating is not damaged. If the coating becomes damaged, discard the electrode. Some of the factors that can increase the likelihood of coating damage are higher power settings, longer activation times, and more arcing.

Do not clean the electrode with a scratch pad or other abrasive object. Scraping with a sharp object or bending may damage the electrode. If damaged, discard the electrode.

Cleaning the Instrument During Use

Important

Wipe the electrode often with moist gauze or other material during surgery.

Precaution

Do not activate the instrument while cleaning the electrode. Injury to operating room personnel may result.

After Surgery

1. After the procedure, turn off the smoke evacuation unit and the electrosurgical generator.
2. Disconnect the pencil, electrode assembly, and the smoke evacuation unit from the generator.

Warning

This product cannot be adequately cleaned and/or sterilized by the user in order to facilitate safe reuse, and is therefore intended for single use. Attempts to clean or sterilize these devices may result in a bio-incompatibility, infection, or product failure risks to the patient.

Precaution

Discard biologically contaminated devices in accordance with your institution's hazardous medical waste and sharps procedures and local regulatory requirements.

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