# **PROXIMATE® RH**

**Rotating Head Skin Stapler**

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
---

**Please read all information carefully.**

Failure to properly follow the instructions may lead to serious surgical consequences.

**Important:** This package insert is designed to provide instructions for use of the PROXIMATE® RH Rotating Head Skin Stapler. It is not a reference to surgical stapling techniques.









PROXIMATE is a trademark of Ethicon Endo-Surgery.

**ETHICON ENDO-SURGERY, LLC**

a Johnson & Johnson company

**Instructions, Instructions, Gebrauchsanweisung, Istruzioni, Instruções, Instrucciones, Gebruiksaanwijzing, Brugsvejledning, Ohje, Οδηγίες, Bruksanvisning, Instrukcja, Utasítások, Návod k použití, Návod, 使用说明**

## Indications

The PROXIMATE RH Rotating Head Skin Stapler has application for routine skin closure in a wide variety of surgical procedures.

## Contraindications

When it is not possible to maintain at least a 5 mm distance from the stapled skin to underlying bones, vessels, or internal organs, the use of staples for skin closure is contraindicated.

## Device Description

The PROXIMATE RH Rotating Head Skin Stapler is a sterile, single patient use instrument designed to deliver rectangular, stainless steel staples for routine wound closure. The head rotates 360° to facilitate placement of the staples.

## MR Conditional

Non-clinical testing has demonstrated the removable staples made of 316L stainless steel alloy in these devices are MR Conditional. A patient with the skin staples can be scanned safely immediately after placement of the staples, under the following conditions:

- Static magnetic field of 3.0 Tesla or less
- Highest spatial magnetic gradient field of 6.5 Tesla/m
- Maximum MR System reported, whole body averaged specific absorption rate (SAR) of 1.3 W/kg for 20 minutes of scanning (per pulse sequence).

## MRI Related Heating

In non-clinical testing, staple lines produced a temperature rise of less than 2°C using the following conditions:

- At 3-Tesla (Magnetom Trio Siemens Medical Solutions MR scanner, software version Numaris/4 syngo MRA30), a maximum MR system-reported whole body averaged SAR of 1.3 W/kg
- 20 minutes of continuous MR scanning (per pulse sequence) using transmit/receive RF body coil.

## Artifact Information

MR image quality may be compromised if the area of interest is in the same area or relatively close to the position of the staple. Therefore, optimization of MR imaging parameters to compensate for the presence of the staples may be necessary.

The worst case signal void size for a staple was:

Pulse Sequence	SE	SE	GRE	GRE
Plane Orientation	Parallel	Perpendicular	Parallel	Perpendicular
Signal Void Size (mm <sup>2</sup> )	1342	1220	1619	1486

## Instructions for Use

Verify compatibility of all instruments and accessories prior to using the instrument.

- 1 Using sterile technique, remove the instrument from the package. To avoid damage, do not flip the instrument into the sterile field.
- 2 **Suggested Eversion Techniques:** With two tissue forceps, pick up each wound edge individually and approximate the edges. (Illustration 1)  
**Or,** with one tissue forcep, pull skin edges together until edges evert. (Illustration 2)  
**Or,** apply tension to either end of the incision, such that the tissue edges begin to approximate themselves. One forcep can be used to ensure that the edges are everted. (Illustration 3)
- 3 Position the instrument over the everted skin edges, aligning the instrument arrow with the center of the incision at a 50° - 60° angle from the skin. (Illustration 4)
- 4 Squeeze the trigger until you touch plastic trigger to plastic handle (plastic to plastic), then release the trigger and move the instrument off the incision. (Illustration 5)

**Alternate Release:** If desired, before releasing the trigger lift up on the instrument. This will help to evert the skin edges, which can then be more easily grasped with the tissue forceps. Release the trigger after the forceps are in place, and repeat the sequence to fire the next staple. (Illustration 6)

**Alternate Technique:** The instrument can also be precocked (partially fired) so that the staple points are visible at the nose of the instrument. This feature, in conjunction with the clear nose and alignment arrow, ensures precise staple placement in the skin. (Illustration 7)

Note: If desired, after the instrument has been precocked one leg of the staple can be hooked onto one side of the tissue. This will aid in drawing the tissue together. This technique may be suitable for attaching skin grafts under moderate tension. (Illustration 8)

#### **Warnings and Precautions**

- Be aware of tissue reactions in patients when exposed to allergens (such as nickel).
- Dispose of all opened instruments whether used or unused.
- This device is packaged and sterilized for single use only. Multiple patient use may compromise the device integrity or create a risk of contamination that, in turn, may result in patient injury or illness.

#### **How Supplied**

The stainless steel material is an alloy of nickel, carbon, iron, molybdenum, silicon, and chromium. The PROXIMATE RH Rotating Head Skin Stapler is supplied sterile and preloaded for single patient use. Discard after use.

#### **Formed Staple Dimensions**

Regular staples have an approximate diameter of 0.53 mm, span of 5.7 mm, and leg length of 3.9 mm. Wide staples have an approximate diameter of 0.58 mm, span of 6.9 mm, and leg length of 3.9 mm.

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©EES, LLC 2014
Ethicon Endo-Surgery (Europe) GmbH