
+ Gripping Surface
Technology

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Please read all information carefully.



Failure to properly follow the instructions may lead to serious surgical consequences, such as leakage or disruption.




Important: This package insert is designed to provide instructions for use of the ECHELON FLEX™ 45 mm Powered Plus Articulating Endoscopic Linear Cutters and the ECHELON FLEX™ 60 mm Powered Plus Articulating Endoscopic Linear Cutters. It is not a reference to surgical techniques.

ETHICON


PART OF THE Johnson & Johnson FAMILY OF COMPANIES

Instructions, Instructions, Gebrauchsanweisung, Istruzioni, Instruções, Instrucciones, Instructies, Brugervejledning, Ohjeet, Οδηγίες, Anvisningar, Instrukcje, Utasítások, Návod k použití, Návod na použitie, Instruksjoner, Talimatlar, Инструкции, Instrukciuni, Instruksi, Hướng dẫn, Kasutusjuhend, Instrukcijas, Instrukcija, Инструкции, Upute, Navodila, Упутства, Упутства, 使用说明



7


8a



8b


9


10


11


12


13


14


15


16

17

18

Indications

The ECHELON, ECHELON ENDOPATH™ and ECHELON FLEX™ families of endoscopic linear cutters and reloads are intended for transection, resection, and/or creation of anastomoses. The instruments have application in multiple open or minimally invasive general, gynecologic, urologic, thoracic, and pediatric surgical procedures. They can be used with staple line or tissue buttressing materials. The instruments may also be used for transection and resection of liver parenchyma (hepatic vasculature and biliary structures), pancreas, kidney and spleen.

Contraindications

- Do not use the instruments on the aorta.
- Do not use the instruments on ischemic or necrotic tissue.
- Do not use any linear cutter on major vessels without making provision for proximal and distal control.
- Tissue thickness should be carefully evaluated before firing any stapler. Refer to the **Reload Product Codes Table** below for tissue compression requirement (closed staple height) for each staple size. If tissue cannot comfortably compress to the closed staple height, or easily compresses to less than the closed staple height, the tissue is contraindicated as it may be too thick or too thin for the selected staple size.
- These instruments are not intended for use when surgical stapling is contraindicated.

Undesirable Side Effects/Residual Risks

Undesirable side effects and risks associated with powered surgical stapling include the potential for bleeding, tissue injury, introduction of non-sterile surfaces or pathogen transfer, inflammatory or unintended tissue reaction, property or environmental damage, electrical shock, and foreign body magnetic resonance incompatibility. Also, unintended harm, extended surgery, or altered surgical approach may result from staple line failures, inability to cut, or damaged devices.

Device Description

The ECHELON FLEX™ Powered Plus Articulating Endoscopic Linear Cutters are sterile, single patient use instruments that simultaneously cut and staple tissue. There are six staggered rows of staples, three on either side of the cut line. The ECHELON FLEX™ 45 Powered Plus instruments have a staple line that is approximately 45 mm long and a cut line that is approximately 42 mm long. The ECHELON FLEX™ 60 Powered Plus instruments have a staple line that is approximately 60 mm long and a cut line that is approximately 57 mm long. The shaft can rotate freely in both directions and an articulation mechanism enables the distal portion of the shaft to pivot to facilitate lateral access to the operative site.

The instruments are packaged with a primary lithium battery pack that must be installed prior to use. There are specific requirements for disposing of the battery pack. Refer to the Battery Pack Disposal section.

The instruments are packaged without a reload and must be loaded prior to use. A staple retaining cap on the reload protects the staple leg points during shipping and transportation. The instruments' lock-out feature is designed to prevent a used or improperly installed reload from being refired or an instrument from being fired without a reload.

Caution: Do not load the instrument more than 12 times for a maximum of 12 firings per instrument. The use of the instrument with staple line reinforcement material may reduce the number of firings.

Instrument Product Codes Table

Instrument Code	Description	Shaft Length
PCEE45A	ECHELON FLEX 45 mm Powered Plus Compact Articulating Endoscopic Linear Cutter	280 mm
PSEE45A	ECHELON FLEX 45 mm Powered Plus Articulating Endoscopic Linear Cutters	340 mm
PLEE45A	ECHELON FLEX 45 mm Powered Plus Long Articulating Endoscopic Linear Cutter	440 mm
PCEE60A	ECHELON FLEX 60 mm Powered Plus Compact Articulating Endoscopic Linear Cutter	280 mm
PSEE60A	ECHELON FLEX 60 mm Powered Plus Articulating Endoscopic Linear Cutter	340 mm
PLEE60A	ECHELON FLEX 60 mm Powered Plus Long Articulating Endoscopic Linear Cutter	440 mm

Reload Product Codes Table

WARNING: The ECHELON FLEX™ 45 mm Powered Plus instruments may only be used with the ECHELON ENDOPATH 45 mm reloads. The ECHELON FLEX™ 60 mm Powered Plus instruments may only be used with the ECHELON ENDOPATH 60 mm reloads.

Reload Code	Open Staple Height	Closed Staple Height	Staple Line Length	Reload Color	Number of Staples	Staple Rows
GST45W	2.6 mm	1.0 mm	45 mm	White	70	6
GST45B	3.6 mm	1.5 mm	45 mm	Blue	70	6
GST45D	3.8 mm	1.8 mm	45 mm	Gold	70	6
GST45G	4.1 mm	2.0 mm	45 mm	Green	70	6
GST45T	4.2 mm	2.3 mm	45 mm	Black	70	6
GST60W	2.6 mm	1.0 mm	60 mm	White	88	6
GST60B	3.6 mm	1.5 mm	60 mm	Blue	88	6
GST60D	3.8 mm	1.8 mm	60 mm	Gold	88	6
GST60G	4.1 mm	2.0 mm	60 mm	Green	88	6
GST60T	4.2 mm	2.3 mm	60 mm	Black	88	6

Clinical Section

Examples of procedures in which linear cutters and staplers may be used include:

General	Gynecologic	Urologic	Thoracic
• Anastomoses	• Vaginal Hysterectomy	• Nephrectomy	• Lung Resection and Biopsy
• Partial Gastropexy		• Cystectomy	• Pulmonary Wedge Resection
• Appendectomy			

Illustrations and Nomenclature (Illustration 1)

- | | | |
|---------------------------------|-----------------------------|---------------------------|
| 1. Closing Trigger | 8. Knife Reverse Switch | 15. Staple Retaining Cap |
| 2. Red Firing Trigger Lock | 9. Rotating Knob | 16. Anvil Jaw |
| 3. Firing Trigger | 10. Articulation Fins | 17. Reload Jaw |
| 4. Anvil Release Button | 11. Reload (Cartridge) | 18. Staple Line |
| 5. Battery Pack | 12. Reload Gripping Surface | 19. Cut Line |
| 6. Battery Pack Release Tab | 13. Reload Alignment Tab | 20. Proximal Black Line |
| 7. Manual Override Access Panel | 14. Reload Alignment Slot | 21. Knife Blade Indicator |

Instructions for Use**Prior to Using Instrument**

- 1 Verify that the reload size matches the instrument size to be used (e.g. use an ECHELON ENDOPATH 45 mm Reload with Gripping Surface Technology reload with an ECHELON FLEX™ 45 Powered Plus instrument).
- 2 Become familiar with the articulating action of the instrument by pulling on the fins of the rotating knob with the index finger and articulating the jaws with the other hand (Illustration 2). As lateral pressure is applied in either direction, the jaws pivot to a maximum angle of 45 degrees, at which point the pivoting force increases to signal the end of the articulating arc. Maintain the lateral pressure while releasing the rotating knob to lock the angle of the shaft. After locking the jaws, release the lateral pressure. The jaws will also lock at approximate 15 and 30 degrees by releasing the rotating knob before the end of the articulation arc. To return the jaws to the straight position, pull on the fins of the rotating knob to release the lock; the jaws will return automatically.
- 3 Verify compatibility of all instruments and accessories prior to using the instrument (refer to **Warnings and Precautions**).

Preparing the Instrument for Use

- 4 Using sterile technique, remove the instrument, battery pack, and reload from their respective packages. To avoid damage, do not flip the instrument, battery pack, or reload into the sterile field.
- 5 Install the battery pack. The battery pack must be installed prior to use. Insert the battery pack by lining up the release tabs on the battery pack with the slots on the rear of the instrument. It can be inserted in either orientation; there is no up or down (Illustration 3). Ensure battery pack is fully inserted into the device. An audible click will be heard when the battery pack is fully inserted.
Caution: The instrument must be used within 12 hours of inserting the battery pack. Refer to the **Battery Pack Disposal** section for battery pack disposal instructions.
- 6 Prior to inserting the reload, ensure the instrument is in the open position (Illustration 1).
- 7 Examine the reload for the presence of a staple retaining cap. If the retaining cap is not in place, discard the reload.
Caution: Tissue thickness should be carefully evaluated prior to using the instrument. Refer to the **Reload Product Codes Table** for proper reload selection.
Caution: Selection of the appropriate staple reload should be based upon the combined thicknesses of both the tissue and the staple line reinforcement material. The use of staple line reinforcement material with the instrument may require an increased force to close and may reduce the number of times the device may be fired. When using staple line reinforcement material, the instructions of the manufacturer of the material should be followed.
- 8 Insert the reload by sliding it against the bottom of the reload jaw until the reload alignment tab stops in the reload alignment slot. Snap the reload securely in place. Remove the staple retaining cap and discard (Illustration 4). The instrument is now loaded and ready for use.
Caution: After removing the staple retaining cap, observe the surface of the new reload. The reload must be replaced with another reload if any colored drivers are visible. (If colored drivers are visible, the reload may not contain staples.)

Using the Instrument

- 9 Close the jaws of the instrument by squeezing the closing trigger until it locks in place (Illustration 5). An audible click indicates that the closing trigger and the jaws are locked. When the jaws of the instrument are closed, the red firing trigger lock and firing trigger are exposed.
Caution: Do not pull the red firing trigger lock or firing trigger at this time. The instrument may be partially or completely fired and will need to be reloaded before using on tissue.
 - 10 Visually inspect the stapler to ensure proper reload seating. Introduce the instrument into the body cavity through a trocar of the appropriate size or through an incision (Illustration 6). When using a trocar, the instrument jaws must be visible past the trocar sleeve before opening the jaws.
Caution: For insertion and removal of articulating instruments, the jaws of the instrument must be straight, parallel to the shaft of the instrument. Failure to have the instrument jaws in the straight position will result in difficult insertion or withdrawal of the instrument and may result in damage to the instrument.
-  **Caution:** When placing the instrument through the trocar or incision, avoid inadvertently pulling the red firing trigger lock and the firing trigger. The instrument may be partially or completely fired and will need to be reloaded before using on tissue. If the instrument is partially fired, remove the instrument and slide the knife reverse switch forward to return the knife to home position (Illustration 7). To open the jaws, squeeze the closing trigger, then simultaneously press the anvil release button on either side of the instrument (Illustration 8a). While pressure is still on the anvil release button, slowly release the closing trigger (Illustration 8b). Remove the reload and insert a new one (see Reloading the Instrument). If the instrument is fired completely, the knife will return to home position automatically. Remove the instrument and insert a new reload (see Reloading the Instrument).

-  a. At any time, if the knife reverse switch does not return the knife to home position and the jaws will not open:
 - First, ensure the battery pack is securely installed and the instrument has power (Illustration 9); then, try the knife reverse switch again (Illustration 7).
 - If the knife still does not return, use the manual override.**Caution:** After the manual override system is used, the instrument is disabled and cannot be used for any subsequent firings. To use the manual override, remove the access panel labeled "Manual Override" on the top of the instrument handle. The manual override lever will be exposed. Move the lever forward and backward until it can no longer be moved (Illustration 10). The knife will now be in the home position. This can be verified by viewing the position of the knife blade indicator on the underside of the reload jaw (Illustration 11). Discard the instrument.
- 11 Once in the cavity, squeeze the closing trigger, then simultaneously press the anvil release button on either side of the instrument (Illustration 8a). While pressure is still on the anvil release button, slowly release the closing trigger to reopen the jaws and return the closing trigger to its original position (Illustration 8b).
- 12 If needed, rotate the jaws by pushing on the fins of the rotating knob with the index finger using a downward or upward pressure (Illustration 12). The instrument shaft will rotate freely in either direction.
- 13 To articulate the jaws inside the body cavity, select an appropriate surface (a body structure, organ, or another instrument) to use as a grounding surface to bend the jaws, ensuring that it stays within the field of view. Rest the lateral side of the jaws that are opposite to the desired direction of articulation against the grounding surface (the jaws must be open in order to articulate the instrument) (Illustration 13).
- WARNING:** Do not attempt to articulate by pressing the front of the jaws against the grounding surface as tissue damage or tissue trauma may occur (Illustration 14). Pull the fins of the rotating knob back with the index finger and apply a sweeping motion towards the side articulation is wanted while gently pushing the instrument handle towards the grounding surface. Maintain the jaws pressed against the grounding surface during this action. Once the desired articulation angle is reached, release the rotating knob to lock the angle (the instrument will only lock at pre-determined angles – 15°, 30° and 45°).
- Caution:** The instrument can only achieve a maximum articulation angle of 45°. When using body structures or organs as a grounding surface, particular attention should be placed to the visual cues and tactile feedback received from the instrument. When the maximum angle is reached, the force will increase indicating the maximum angle has been reached. Avoid applying excessive pressure to the tissue as tissue damage or tissue trauma may occur.
- 14 Position the instrument around the tissue to be stapled.
- Caution:** Ensure that the tissue lies flat and is positioned properly between the jaws. Any "bunching" of tissue along the reload, particularly in the crotch of the jaws, may result in an incomplete staple line. The distal black lines on the anvil and reload channel designate the ends of the staple line. The line on the reload channel that reads "cut" references the cut line on the device (Illustration 1).
- Caution:** When positioning the stapler on the application site, ensure that no obstructions such as clips, stents, guide wires, etc. are within the instrument jaws. Firing over an obstruction may result in incomplete cutting action, improperly formed staples, and/or inability to open the instrument jaws.
- 15 After positioning the instrument jaws, close the jaws by squeezing the closing trigger until it locks (Illustration 5). An audible click indicates that the closing trigger and the jaws are locked. When the jaws of the instrument are closed, the red firing trigger lock and firing trigger are exposed. Holding the jaws in place for 15 seconds after closing and prior to firing may result in better compression and staple formation.
-  **Caution:** Ensure tissue has not extended (extruded) proximal to the proximal black line on the instrument (Illustration 1). Tissue forced into the instrument proximal to the black line may be transected without staples.
- Caution:** If the closing trigger is difficult to lock, **reposition the instrument and take a smaller amount of tissue**. Ensure that the proper reload selection has been made. (Refer to the **Reload Product Codes Table**.)
- Caution:** If the clamping mechanism becomes inoperative and the jaws do not clamp on tissue, do not fire the instrument. Remove and do not continue to use instrument.
- Caution:** The use of staple line buttressing materials with the instrument may require an increased force to close. When using staple line reinforcement material, the instructions of the manufacturer of the material should be followed.
- 16 Pull back the red firing-trigger lock to enable the firing trigger to be pulled (Illustration 15).
- 17 Fire the instrument by pulling the firing trigger; the motor will activate audibly (Illustration 16). Continue to depress the trigger until the motor stops (audible feedback). The motor will slow down when device encounters more force. Alternatively, for enhanced tissue compression, the user may pulse the firing trigger throughout the firing cycle until the transection is completed.
- Caution:** Since the motor may stop if it stalls, it is important to do a visual check to ensure that the knife blade indicator, on the underside of the reload jaw, has reached the end of the transection.
-  **Caution:** Attempting to force the device to complete the firing stroke with too much tissue between the jaws, or with dense/thick tissue between the jaws, may result in motor stall and the knife will stop. If this happens, release the firing trigger, slide the knife reverse switch forward, and remove and reload the instrument (Illustration 7). Then, position the instrument around a smaller tissue section or use a more appropriate reload (refer to the **Reload Product Codes Table**).
- 18 To complete the firing sequence, release the firing trigger to activate the motor and automatically return the knife to home position where the motor will stop. In this position, the instrument is locked out until the jaws are opened and re-closed.
 - If it is necessary to interrupt the firing sequence or it is interrupted inadvertently by releasing the trigger during the firing sequence, pull the firing trigger again to continue. The status of the transection can be determined by observing the knife blade indicator on the underside of the reload jaw at the beginning and end of the firing (Illustration 1). Once the knife reaches the end of the firing, releasing the trigger automatically returns the knife to home position.

 **Caution:** If the instrument locks out, the motor will stop. Release the firing trigger and slide the knife reverse switch forward to return the knife to the home position (Illustration 7). In this position, the instrument should be removed, opened, and reloaded in order to continue. To open the jaws, squeeze the closing trigger, then simultaneously press the anvil release button on either side of the instrument (Illustration 8a). While pressure is still on the anvil release button, slowly release the closing trigger (Illustration 8b). Follow the instructions for Reloading the Instrument.

 b. At any time, if the knife reverse switch does not return the knife to home position and the jaws will not open:

 - First, ensure the battery pack is securely installed and the instrument has power (Illustration 9). Then, try the knife reverse switch again (Illustration 7).
 - If the knife still does not return, use the manual override.**Caution:** After the manual override system is used, the instrument is disabled and cannot be used for any subsequent firings. To use the manual override, remove the access panel labeled "Manual Override" on the

top of the instrument handle. The manual override lever will be exposed. Move the lever forward and backward until it can no longer be moved (Illustration 10). The knife will now be in the home position. This can be verified by viewing the position of the knife blade indicator on the underside of the reload jaw (Illustration 11). Discard the instrument.

Caution: Incomplete firing may result in malformed staples, incomplete cut line, bleeding, and/or difficulty removing the device.

Caution: The use of staple line reinforcement material with the instrument may reduce the number of times the device may be fired. When using staple line reinforcement material, the instructions of the manufacturer of the material should be followed.

Caution: Crossing of staple lines may shorten the life of the instrument.

Caution: If the firing mechanism becomes inoperative, **do not** continue to use the instrument.

- 19 To open the jaws, squeeze the closing trigger, then simultaneously press the anvil release button on either side of the instrument. While pressure is still on the anvil release button, slowly release the closing trigger. (Illustration 8).

Caution: If the jaws do not automatically open after the anvil release button is pressed, first ensure that the knife is in the home position. The position of the knife can be determined by observing the knife blade indicator under the reload jaw (Illustration 11). If the knife blade indicator is not in the home position or the position of the knife cannot be determined, slide the knife return switch to activate the motor and return the knife to home position (Illustration 7). Try opening the jaws again using the anvil release button. If the jaws do not open at this point, then gently pull the closing trigger (1) upward (away from the handle) until both firing and closing triggers return to their original positions.
- 20 Gently pull the instrument away from the transected tissue and ensure it is released from the jaws.

Caution: Examine the staple lines for pneumostasis/hemostasis and proper staple closure. Minor bleeding can be controlled with manual sutures or other appropriate techniques.
- 21 Before removing articulating instruments, move the jaws away from any obstruction inside the body cavity while keeping the jaws open and within the field of view, and pull on the fins of the rotating knob. The jaws will return to the straight position automatically.

Caution: For insertion and removal of articulating instruments, the jaws of the instrument must be straight, parallel to the shaft of the instrument. Failure to have the instrument jaws in the straight position will result in difficult insertion or withdrawal of the instrument and may result in damage to the instrument.
- 22 To remove the instrument from the cavity, squeeze the closing trigger until it locks, closing the jaws (Illustration 5). Completely withdraw the instrument in the closed position.

Reloading the Instrument

23 Using sterile technique, remove the reload from the package. To avoid damage, do not flip the reload into the sterile field.

24 Prior to reloading, ensure the instrument is in the open position (Illustration 1).

25 Push upward (toward the anvil) to unsnap the reload from the reload jaw. Discard the used reload (Illustration 17).

 **WARNING:** Prior to reloading the instrument, hold the instrument in a vertical position, with anvil and reload jaw completely submerged in sterile solution. Swish vigorously and then wipe the inside and outside surfaces of the anvil and reload jaw to clean any unused staples from the instrument. Do not use the instrument until it has been visually inspected to confirm there are no staples on the anvil and reload jaw.

26 Examine the new reload for the presence of a staple retaining cap. If the retaining cap is not in place, discard the reload.

Caution: Tissue thickness should be carefully evaluated prior to using the instrument. Refer to the **Reload Product Codes Table** for proper reload selection.

Caution: Selection of the appropriate staple reload should be based upon the combined thicknesses of both the tissue and the staple line reinforcement material. The use of staple line reinforcement material with the instrument may require an increased force to close and may reduce the number of times the device may be fired. When using staple line reinforcement material, the instructions of the manufacturer of the material should be followed.

27 Insert the new reload by sliding it against the bottom of the reload jaw until the reload alignment tab stops in the reload alignment slot. Snap the reload securely in place. Remove the staple retaining cap and discard. The instrument is now reloaded and ready for use (Illustration 4).

Caution: After removing the staple retaining cap, observe the surface of the new reload. The reload must be replaced with another reload if any colored drivers are visible. (If colored drivers are visible, the reload may not contain staples.)

Battery Pack Disposal

The battery pack contains a built-in battery drain and must be **installed** into the instrument to trigger the automatic battery drain feature. It does not need to remain in the instrument—it will continue to drain after it has been removed from the instrument. Once removed, depending on your local regulations, it can be put directly into the battery recycling bin or normal waste stream. Batteries used in this device contain 1,2-dimethoxyethane.

Prior to Installation into the Instrument

If battery pack needs to be disposed of prior to installation into instrument (e.g. product is beyond expiration date indicated on the package, battery pack is dropped), first install the battery pack into the instrument, then remove, to trigger built-in battery drain.

After Instrument Use

The battery pack must be removed from the instrument prior to disposal.

If battery pack requires decontamination prior to disposal, follow hospital protocol or the **Battery Pack Cleaning and Disinfection** instructions below.

To remove the battery pack, squeeze the release tabs and pull the battery pack straight back (Illustration 18). It is not necessary to disassemble the battery pack.

Battery Pack Cleaning and Disinfection

WARNING: Do not use an autoclave, ethylene oxide, or radiation to sterilize or disinfect the battery pack.

Manual Cleaning

- 1 Remove the battery pack from the instrument before cleaning. Battery pack should not be submerged in water or cleaning solutions.
- 2 Clean the battery pack surfaces with a neutral pH detergent or neutral pH enzymatic detergent, prepared according to the manufacturer's instructions.
- 3 Use soft bristle brush to manually clean the battery pack with the cleaning solution.
- 4 Ensure areas containing crevices are scrubbed thoroughly.
- 5 Wipe off detergent thoroughly with lukewarm tap water.
- 6 Perform visual inspection to determine if debris is removed.
- 7 Repeat cleaning as necessary to obtain a visually clean battery pack.

Chemical Disinfection

Disinfectants should be prepared and used according to the manufacturer's recommendations. It is recommended that the chemical disinfectant be wiped off with tap water.

- Simple Green D Pro 3®
- 10% Bleach (sodium hypochlorite solution)
- 70% Isopropyl alcohol

Standard Conventions Used

The Use of Caution, Warning, and Note Statements

Information relative to the completion of a task in a safe and thorough manner will be supplied in the form of a Caution or a Warning statement. These statements are found throughout the documentation.

These statements should be read before continuing to the next step in a procedure.

WARNING: A Warning statement indicates an operating or maintenance procedure, practice, or condition that, if not strictly observed, could result in personal injury or loss of life.

Caution: A Caution statement indicates an operating or maintenance procedure, practice, or condition that, if not strictly observed, could result in damage to or destruction of the equipment.

Warnings and Precautions

- Minimally invasive procedures should be performed only by persons having adequate training and familiarity with minimally invasive techniques. Consult medical literature relative to techniques, complications, and hazards prior to performance of any minimally invasive procedure. Failure to properly follow the instructions may lead to serious surgical consequences, such as leakage or disruption.
- Minimally invasive instruments may vary in diameter from manufacturer to manufacturer. When minimally invasive instruments and accessories from different manufacturers are employed together in a procedure, verify compatibility prior to initiation of the procedure.
- When using other technologies (e.g., electrocautery) in the procedure, observe the precautions suggested by the original equipment manufacturer to avoid the hazards associated with their use.
- Pre-operative radiotherapy may result in changes to tissue. These changes may, for example, cause the tissue thickness to exceed the indicated range for the selected staple. Careful consideration should be given to any pre-surgical treatment the patient may have undergone, which may require alterations to surgical technique or alternative surgical procedures.
- Do not use the ECHELON FLEX™ Powered Plus Articulating Endoscopic Linear Cutters after the expiration date or if the package is damaged. This may render the equipment inoperable or non-sterile.
- Do not load the instrument more than 12 times for a maximum of 12 firings per instrument. The use of the instrument with staple line reinforcement material may reduce the number of firings.
- The ECHELON FLEX™ 45 mm Powered Plus instruments may only be used with ECHELON ENDOPATH 45 mm reloads. The ECHELON FLEX™ 60 mm Powered Plus instruments may only be used with the ECHELON ENDOPATH 60 mm reloads.
- The instrument must be used within 12 hours of inserting the battery pack. Refer to the **Battery Pack Disposal** section for battery pack disposal instructions.
- Tissue thickness should be carefully evaluated prior to using the instrument. Refer to the **Reload Product Codes Table** for proper reload selection.
- Selection of the appropriate staple reload should be based upon the combined thicknesses of both the tissue and the staple line reinforcement material. The use of staple line reinforcement material with the instrument may require an increased force to close and may reduce the number of times the device may be fired. When using staple line reinforcement material, the instructions of the manufacturer of the material should be followed.
- After removing the staple retaining cap, observe the surface of the new reload. The reload must be replaced with another reload if any colored drivers are visible. (If colored drivers are visible, the reload may not contain staples.)
- Do not pull the red firing trigger lock or firing trigger at this time. The instrument may be partially or completely fired and will need to be reloaded before using on tissue.
- For insertion and removal of articulating instruments, the jaws of the instrument must be straight, parallel to the shaft of the instrument. Failure to have the instrument jaws in the straight position will result in difficult insertion or withdrawal of the instrument and may result in damage to the instrument.
- When placing the instrument through the trocar or incision, avoid inadvertently pulling the red firing trigger lock and the firing trigger. The instrument may be partially or completely fired and will need to be reloaded before using on tissue.

- Do not attempt to articulate by pressing the front of the jaws against the grounding surface as tissue damage or tissue trauma may occur.
- The instrument can only achieve a maximum articulation angle of 45°. When using body structures or organs as a grounding surface, particular attention should be placed to the visual cues and tactile feedback received from the instrument. When the maximum angle is reached, the force will increase indicating the maximum angle has been reached. Avoid applying excessive pressure to the tissue as tissue damage or tissue trauma may occur.
- Ensure that the tissue lies flat and is positioned properly between the jaws. Any “bunching” of tissue along the reload, particularly in the crotch of the jaws, may result in an incomplete staple line.
- When positioning the stapler on the application site, ensure that no obstructions such as clips, stents, guide wires, etc. are within the instrument jaws. Firing over an obstruction may result in incomplete cutting action, improperly formed staples, and/or inability to open the instrument jaws.
- Ensure tissue has not extended (extruded) proximal to the proximal black line on the instrument. Tissue forced into the instrument proximal to the black line may be transected without staples.
- If the closing trigger is difficult to lock, **reposition the instrument and take a smaller amount of tissue**. Ensure that the proper reload selection has been made. (Refer to the **Reload Product Codes Table**.)
- If the clamping mechanism becomes inoperative and the jaws do not clamp on tissue, do not fire the instrument. Remove and do not continue to use instrument.
- The use of staple line buttressing materials with the instrument may require an increased force to close. When using staple line reinforcement material, the instructions of the manufacturer of the material should be followed.
- Attempting to force the trigger to complete the firing stroke with too much tissue between the jaws, or with dense/thick tissue between the jaws, may result in motor stall and the knife will stop. If this happens, release the firing trigger, slide the knife reverse switch forward, and remove and reload the instrument (Illustration 7). Then, position the instrument around a smaller tissue section or use a more appropriate reload (refer to the **Reload Product Codes Table**).
- Since the motor may stop if it stalls, it is important to do a visual check to ensure that the knife blade indicator, on the underside of the reload jaw, has reached the end of the transection.
- If the instrument locks out, the motor will stop. Release the firing trigger and slide the knife reverse switch forward to return the knife to the home position. In this position, the instrument should be removed, opened, and reloaded in order to continue.
- After the manual override system is used, the instrument is disabled and cannot be used for any subsequent firings. To use the manual override, remove the access panel labeled “Manual Override” on the top of the instrument handle. The manual override lever will be exposed. Move the lever forward and backward until it can no longer be moved (Illustration 10). The knife will now be in the home position. This can be verified by viewing the position of the knife blade indicator on the underside of the reload jaw (Illustration 11). Discard the instrument.
- Selection of the appropriate staple reload should be based upon the combined thicknesses of both the tissue and the staple line reinforcement material. The use of staple line reinforcement material with the instrument may reduce the number of times the device may be fired. When using staple line reinforcement material, the instructions of the manufacturer of the material should be followed.
- Incomplete firing may result in malformed staples, incomplete cut line, bleeding, and/or difficulty removing the device.
- Crossing of staple lines may shorten the life of the instrument.
- If the firing mechanism becomes inoperative, **do not** continue to use the instrument.
- If the jaws do not automatically open after the anvil release button is pressed, first ensure that the knife is in the home position. The position of the knife can be determined by observing the knife blade indicator under the reload jaw. If the knife blade indicator is not in the home position or the position of the knife cannot be determined, slide the knife return switch to activate the motor and return the knife to home position. Try opening the jaws again using the anvil release button. If the jaws do not open at this point, then gently pull the closing trigger upward (away from the handle) until both firing and closing triggers return to their original positions.
- Examine the staple lines for pneumostasis/hemostasis and proper staple closure. Minor bleeding can be controlled with manual sutures or other appropriate techniques.
- Prior to reloading the instrument, hold the instrument in a vertical position, with anvil and reload jaw completely submerged in sterile solution. Swish vigorously and then wipe the inside and outside surfaces of the anvil and reload jaw to clean any unused staples from the instrument. Do not use the instrument until it has been visually inspected to confirm there are no staples on the anvil and reload jaw.
- Before removing the instrument, be sure tissue is cleared from the jaws and then close the jaws.
- When selecting the reload, careful consideration should be given to existing pathologic conditions as well as any pre-surgical treatment, such as radiotherapy, that the patient may have undergone. Certain conditions or preoperative treatments may cause change in tissue thickness that would exceed the indicated range of tissue thickness for the standard choice of reload.
- When dividing major vascular structures, be sure to adhere to the basic surgical principle of proximal and distal control.
- Do not use an autoclave, ethylene oxide, or radiation to sterilize or disinfect the battery pack.
- The battery may present a fire hazard if mistreated. Do not disassemble, heat above 100 °C, autoclave, crush, puncture, short external contacts or recharge.
- Use of any other type of battery other than the battery supplied with the device may result in increased EMISSIONS or decreased IMMUNITY of the ECHELON FLEX™ Powered Plus Articulating Endoscopic Linear Cutters.
- Portable radio and mobile RF communications equipment can affect medical devices. Follow the directions of the following tables when using the ECHELON FLEX™ Powered Plus Articulating Endoscopic Linear Cutters.
- Avoid use of the ECHELON FLEX™ Powered Plus Articulating Endoscopic Linear Cutters adjacent to or stacked with other equipment. If adjacent or stacked use is necessary, monitor the ECHELON FLEX™ Powered Plus Articulating Endoscopic Linear Cutters and the other equipment to ensure normal operation.
- Do not modify this equipment without authorization from the manufacturer.
- Instruments or devices which come into contact with bodily fluids may require special disposal handling to prevent biological contamination.
- Dispose of all opened devices whether used or unused. This device is packaged and sterilized for single use only.
- Reuse and improper reprocessing or resterilization may compromise the structural integrity of the device and/or lead to device failure that in turn may result in patient injury, illness, or death.
- Reuse and improper reprocessing or resterilization of single use devices may create a risk of contamination and/or cause infection or cross-infection, including but not limited to, the transmission of infectious diseases. Contamination may lead to injury, illness, or death.
- Federal (USA) law restricts this device to sale by or on the order of a licensed healthcare practitioner.

Specifications

The ECHELON FLEX™ Powered Plus Articulating Endoscopic Linear Cutters have a power rating of 40 Watts.

The ECHELON FLEX™ Powered Plus Articulating Endoscopic Linear Cutters are resistant to water ingress and are classified per IEC 60601-1 as IPX0.

The ECHELON FLEX™ Powered Plus Articulating Endoscopic Linear Cutters need special precautions regarding EMC and need to be installed and put into service according to the EMC information provided in this document. Portable and mobile RF communications equipment can affect Medical Electrical Equipment.

WARNING: Use of any other type of battery other than the battery supplied with the device may result in increased EMISSIONS or decreased IMMUNITY of the ECHELON FLEX™ Powered Plus Articulating Endoscopic Linear Cutters.

WARNING: Portable radio and mobile RF communications equipment can affect Medical Electrical Equipment. Follow the directions of the following tables when using the ECHELON FLEX™ Powered Plus Articulating Endoscopic Linear Cutters.

WARNING: Avoid use of the ECHELON FLEX™ Powered Plus Articulating Endoscopic Linear Cutters adjacent to or stacked with other equipment. If adjacent or stacked use is necessary, monitor the ECHELON FLEX™ Powered Plus Articulating Endoscopic Linear Cutters and the other equipment to ensure normal operation.

Table of Toxic and Hazardous Substances / Elements and their Content

As required by China's Management Methods for Controlling Pollution

By Electronic Information Products

Product Codes: ECHELON FLEX™ POWERED ENDOPATH series

	Toxic and Hazardous Substances and Elements					
Part Name	Lead (Pb)	Mercury (Hg)	Cadmium (Cd)	Hexavalent Chromium (Cr(VI))	Polybrominated biphenyls (PBB)	Polybrominated diphenyl ethers (PBDE)
Handle	X	O	O	O	O	O
Shaft	O	O	O	O	O	O

O: indicates that this toxic or hazardous substance contained in all of the homogenous materials for this part is below the limit requirement in SJ/T11363-2006.

X: indicates that this toxic or hazardous substance contained in at least one of the homogeneous materials used for this part is above the limit requirement in SJ/T11363-2006.

All parts named in this table are in compliance with the European union's RoHS legislation "Directive 2011/65/EU of the European Parliament and of the Council of 8 June 2011 on the restriction of the use of hazardous substances in electrical and electronic equipment."

The referenced Environmental Protection Use Period Marking was determined according to normal operating conditions of the product such as temperature and humidity.

Electromagnetic Emissions Guidance

Guidance and manufacturer's declaration—electromagnetic emissions		
The ECHELON FLEX™ Powered Plus Articulating Endoscopic Linear Cutters are intended for use in the electromagnetic environment specified below. The customer or the user of the ECHELON FLEX™ Powered Plus Articulating Endoscopic Linear Cutters should ensure that they are used in such an environment.		
Emissions test	Compliance	Electromagnetic environment - guidance
Harmonic emissions IEC 61000-3-2	Not applicable	The ECHELON FLEX™ Powered Plus Articulating Endoscopic Linear Cutter is suitable for use in all establishments other than domestic and those directly connected to the public low-voltage power supply network that supplies buildings used for domestic purposes.
Voltage fluctuations/ flicker emissions IEC 61000-3-3	Not applicable	
RF emissions CISPR 14-1	Complies	The ECHELON FLEX™ Powered Plus Articulating Endoscopic Linear Cutters are not suitable for interconnection with other equipment.

Electromagnetic Immunity Guidance

For electromagnetic immunity, essential performance is: no knife movement without activation of the firing trigger.

Guidance and manufacturer's declaration—electromagnetic immunity

The ECHELON FLEX™ Powered Plus Articulating Endoscopic Linear Cutters are intended for use in the electromagnetic environment specified below. The customer or the user of the ECHELON FLEX™ Powered Plus Articulating Endoscopic Linear Cutters should ensure that they are used in such an environment.

IMMUNITY test	IEC 60601 test level	Compliance level	Electromagnetic environment — guidance
Electrostatic discharge (ESD) IEC 61000-4-2	± 6 kV contact ± 8 kV air	± 6 kV contact ± 8 kV air	Floors should be wood, concrete or ceramic tile. If floors are covered with synthetic material, the relative humidity should be at least 30 %.
Electrical fast transient/burst IEC 61000-4-4	± 2 kV for power supply lines ± 1 kV for input/output lines	Not applicable	
Surge IEC 61000-4-5	± 1 kV line(s) to line(s) ± 2 kV line(s) to earth	Not applicable	
Voltage dips, short interruptions and voltage variations on power supply input lines IEC 61000-4-11	<5 % U_T (>95 % dip in U_T) for 0.5 cycles 40 % U_T (60 % dip in U_T) for 5 cycles 70 % U_T (30 % dip in U_T) for 25 cycles <5 % U_T (>95 % dip in U_T) for 5 s	Not applicable	Not applicable
Power frequency (50/60 Hz) magnetic field IEC 61000-4-8	3 A/m	3 A/m	Power frequency magnetic fields should be at levels characteristic of a typical location in a typical commercial or hospital environment.

NOTE U_T is the a.c. mains voltage prior to application of the test level.

Guidance and manufacturer's declaration—electromagnetic immunity

The ECHELON FLEX™ Powered Plus Articulating Endoscopic Linear Cutters are intended for use in the electromagnetic environment specified below. The customer or the user of the ECHELON FLEX™ Powered Plus Articulating Endoscopic Linear Cutters should ensure that they are used in such an environment.

IMMUNITY test	IEC 60601 TEST LEVEL	Compliance level	Electromagnetic environment – guidance
			<p>Portable and mobile RF communications equipment should be used no closer to any part of the Powered Plus Articulating Endoscopic Linear Cutter, including cables, than the recommended separation distance calculated from the equation applicable to the frequency of the transmitter.</p> <p>Recommended separation distance</p> <p>Conducted RF IEC 61000-4-6 3 Vrms 150 kHz to 80 MHz Not applicable</p> <p>Radiated RF IEC 61000-4-3 3 V/m 80 MHz to 2.5 GHz 3 V/m</p> <p>Not Applicable</p> $d = 1.2\sqrt{P} \quad 80 \text{ MHz to } 800 \text{ MHz}$ $d = 2.3\sqrt{P} \quad 800 \text{ MHz to } 2.5 \text{ GHz}$ <p>where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer and d is the recommended separation distance in meters (m).</p> <p>Field strengths from fixed RF transmitters, as determined by an electromagnetic site survey,^a should be less than the compliance level in each frequency range.</p> <p>Interference may occur in the vicinity of equipment marked with the following symbol:</p> 
NOTE 1 At 80 MHz and 800 MHz, the higher frequency range applies.			

NOTE 2 These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.

^aField strengths from fixed transmitters, such as base stations for radio (cellular/cordless) telephones and land mobile radios, amateur radio, AM and FM radio broadcast and TV broadcast cannot be predicted theoretically with accuracy. To assess the electromagnetic environment due to fixed RF transmitters, an electromagnetic site survey should be considered. If the measured field strength in the location in which the ECHELON FLEX™ Powered Plus Articulating Endoscopic Linear Cutters are used exceeds the applicable RF compliance level above, the ECHELON FLEX™ Powered Plus Articulating Endoscopic Linear Cutters should be observed to verify normal operation. If abnormal performance is observed, additional measures may be necessary, such as re-orienting or relocating the ECHELON FLEX™ Powered Plus Articulating Endoscopic Linear Cutters

Recommended separation distances between portable and mobile RF communications equipment and the ECHELON FLEX™ Powered Plus Articulating Endoscopic Linear Cutters			
Rated maximum output power of transmitter W	Separation distance according to frequency of transmitter m		
	150 kHz to 80 MHz	80 MHz to 800 MHz $d = 1.2\sqrt{P}$	800 MHz to 2.5 GHz $d = 2.3\sqrt{P}$
0.01	Not applicable	0.12	0.23
0.1		0.38	0.73
1		1.2	2.3
10		3.8	7.3
100		12	23

For transmitters rated at a maximum output power not listed above, the recommended separation distance d in meters (m) can be estimated using the equation applicable to the frequency of the transmitter, where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer.

NOTE 1 At 80 MHz and 800 MHz, the separation distance for the higher frequency range applies.

NOTE 2 These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.

Transport and Storage Conditions

Temperature: -20 °C — 25 °C

Relative Humidity: 10% — 60%

Pressure: 500 — 1060 hPa

Operating Environment Conditions

Temperature: 10 °C — 40 °C

Relative Humidity: 30% — 75%

Pressure: 800 — 1060 hPa

How Supplied

The ECHELON FLEX™ Powered Plus Articulating Endoscopic Linear Cutters are supplied sterile for single patient use. Discard after use.

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©Ethicon Endo-Surgery, Inc. 2020
Ethicon Endo-Surgery (Europe) GmbH
ETHICON ENDO-SURGERY, LLC