TME TEMPORARY MYOCARDIAL HEARTWIRES

Instructions for Use

**Temporary Myocardial Heartwires
Unipolar, Bipolar & Quadripolar**

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








en-Quantity
de-Menge
es-Cantidad
fr-Quantité
it-Quantità
pt-Quantidade

English | Table of Contents

1.	Intended Use.....	7
2.	Device Description.....	7
3.	Cautions.....	7
4.	Contraindications.....	7
5.	Adverse Events.....	7
6.	Potential Complications.....	7
7.	Precautions.....	8
8.	Directions for Use.....	8
9.	Handling and Storage.....	9
10.	Limited Warranty and Disclaimer.....	9

Copy - For Information Only

1. Intended Use:

Unipolar, bipolar, and quadripolar temporary myocardial heartwires are intended for temporary atrial and ventricular pacing and sensing during or after cardiac surgery in conjunction with an external pulse generator.

2. Device Description:

Temporary myocardial heartwire, model TME, consists of insulated wire with a curved or straight needle on one side and a heart needle or V-hooks on the other side. The fixation is either tine, zigzag or suture-on. Unipolar=1 pin adaptor set. Bipolar=2 pin adaptor set. Quadripolar= 4 pin adaptor set.

3. Cautions:

- Device is supplied sterile. Do not use if package has been previously opened or damaged.
- Prior to use, read all package inserts, warnings, precautions, and instructions. Failure to do so may result in severe patient injury or death.
- Procedure must be performed by trained medical personnel well versed in anatomical landmarks, safe technique, and potential complications.
- The removal or extraction of the heartwire requires utmost care and should be performed by an experienced physician only.
- Line powered equipment used in the vicinity of the patient must be properly grounded.
- Avoid the possibility of unintentional contact between the heartwires, including extension cable, and any equipment used as well as any conductive surface contact.
- Leaving the heartwire implanted for longer than 7 days may result in difficulty or inability to extract the lead and/or bleeding.
- The device is designed for single use only.
- Do not resterilize or reuse. Do not alter the device in any way.
- If you have any questions, please call Oscor Inc. at (727) 937-2511.

4. Contraindications:

Do not use when a pacemaker has been implanted.

5. Adverse Events:

Heartwire related problems include, but are not limited to:

Event Possible	Result
Fracture.....	Periodic or continued loss of sensing and/or pacing
Dislodgement.....	Periodic or continued loss of sensing and/or pacing. May require surgical repositioning or removal.
Cardiac Perforation.....	Periodic or continued loss of sensing and/or pacing, diaphragmatic muscle and phrenic nerve stimulation; cardiac tamponade
Myocardial Irritability.....	PVC; ventricular tachycardia; fibrillation
Threshold Elevation	Loss of capture and/or loss of sensing
Infection.....	Heartwire may require surgical removal

6. Potential Complications:

The possibility of dislodgement, lead fractures, threshold elevation and sensing inconsistency exists.

Copy - For Information Only

7. Precautions:

Defibrillation and life support equipment and supplies should be placed on standby and be kept immediately available during lead manipulation, pulse generator connection and removal.

8. Directions for Use:

a. (1) Unipolar Heartwire or Quadripolar (4x Unipolar) Heartwire- The distal end of the unipolar temporary heartwire is embedded in the desired portion of the myocardium (atrial or ventricular) using the small curved needle. It is important that the electrode (zig-zag wire or the exposed wire just before the tines) be implanted in the myocardium, i.e., the electrode should not be visible. See Figures a and b.

Fig. a
Unipolar Heartwire-
Zig-Zag



Fig. b
Unipolar Heartwire-
Tines

For Tined Heartwires

Pull the heartwire through the myocardium until the tines are exposed. Pull back on the heartwire body until the tines are stabilized on top of the epicardium.

a. (2) Bipolar Heartwire or Quadripolar (2x Bipolar) Heartwire- The distal bipolar end of the temporary pacing wire is embedded in the desired portion of the myocardium (atrial or ventricular) using the small curved needle. Both electrodes should not be visible. The preferred technique employs two passes of the needle resulting in the insulation between electrodes being visible. See Figures c and d.

Fig. c
Bipolar Heartwire-
Zig-Zag



Fig. d
Bipolar Heartwire-
Tines

For Tined Heartwires

Pull the heartwire through the myocardium until the tines are exposed. Pull back on the heartwire body until the tines are stabilized on top of the epicardium.

Note: The Physician may elect to vary the fixation procedure in accordance with judgment and experience.

Caution: The temporary heartwires must be placed in such a way that moderate traction will remove it when appropriate.

a. (3) Suture-on Quadripolar (4 x Unipolar) Heartwire

Tie the suture inside the V-hook of the heartwire. See Figure e.

Caution: Only use permanent surgical material to always ensure safe fixation of the wire.



Fig. e

Copy - For Information Only

Skip to step c. for further instructions.

- b. Cut off small curved needle and the protruding wire (leaving approximately 5 mm to 10 mm exposed wire) at epicardial surface.
- c. Pull large needle through the thorax. Cut off needle leaving a minimum of 15 mm (9/16") of bare wire. Connect the enclosed 2 mm pin/adaptor. See Steps d.1 through d.3.
- d. Join each proximal wire to connecting pin as follows:


d.1 Expose pin hole and insert bare wire through hole up to insulation. Bend wires back on pin.

d.2 Push pin cover forward onto pin to completely cover bare wires..

d.3 To cover connecting pin when not in use, push cover over pin. To remove, use a slight twisting motion, taking care not to break the tab.

Caution: Bare connectors and bare wires constitute a direct path to the heart. Small currents (approximately 10 μ A) may be sufficient to cause fibrillation. Only trained personnel should handle heartwires and connecting pins.

- e. Use extension cables to connect proximal ends of the heartwires to external device. The quadripolar heartwire is 220 cm and does not require the use of an extension cable.
- f. Recommended threshold range is from: 0.5 V to 5.0 V.
- g. ECG measurements should be done only with ECG recorders having a floating input. If this is not followed, there is the possibility of fibrillation.

9. Handling and Storage:

Avoid subjecting the device to unusual stresses. The device should be stored at temperatures between 5° C and 30° C.

10. Limited Warranty and Disclaimer:

Oscor Inc. hereby warrants that if any Oscor Inc. product fails to perform within normal tolerances due to a defect in materials or workmanship, Oscor Inc. will provide, at no charge, a replacement of Oscor Inc. product. This limited warranty applies only if each of the following conditions are met:

- A. The product was designed, manufactured and distributed by Oscor Inc.
- B. The product has not been mishandled, reprocessed, or altered in any way.
- C. The product was used before the "USE BEFORE" date marked on the packaging of the product.
- D. The failed product must be returned to Oscor Inc. and becomes the property of Oscor Inc.

No representation or warranty is made that an Oscor Inc. product will not fail. Oscor Inc. disclaims responsibility for any medical complications, including death, resulting from the use of its products. Except as expressly provided by this limited warranty, Oscor Inc. is not

Copy - For Information Only

responsible for any direct, incidental, or consequential damages based on any defect, failure, or malfunction of its products whether the claim is based on warranty, contract, tort, or otherwise. Some states do not allow the exclusion or limitation of consequential damages and so the above limitation of exclusion may not apply to you.

Except as expressly provided by the limited warranty, Oscor Inc. makes no warranty, express or implied, including but not limited to, any implied warranty of merchantability, or fitness for a particular purpose.

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Palm Harbor, FL 34683 U.S.A.