

REF 473105
PRODUCT CODE ET125B

Smallbore Extension Set
with Female Luer Lock Connector and
Spin-Lock® Connector

13 **LENGTH** 13 in. (33 cm) (approx.)
PRIMING VOLUME 0.4 mL (approx.)
GTIN: 04046964185090

STERILE Contents sterile in the only
NONPYROGENIC unopened, undamaged package.
Not made with natural rubber latex or DEHP.
P-2514-7 REV. 1Q15

② For single use only. Do not resterilize.

PRECAUTION: Not intended for use with power injectors.

Instructions for Use:
Use Aseptic Technique.
• Replace set per CDC guidelines and/or institutional protocol.

1. Firmly attach connectors to desired connections.
2. Prime extension set to ensure air is expelled.

Pictorial is for reference only.

B. BRAUN
B. Braun Medical Inc.
Bethlehem, PA 18018-3524 USA
Assembled in Dominican Republic.



2026-03-31
LOT 0061778633