

# CATH•SECURE® EXTENDED TAB

**Product: CATH•SECURE® EXTENDED TAB**  
**Item No. 5445-3**



**Item No.** 5445-3

**Item Name:** CATH•SECURE® Extended Tab

**Pack:**

50 each per box

8 box per case

400 each per case

**Description:**

Urological Foley, Dialysis, PICC, Multiple Purpose Securement Device

**Measurements:** 2.75" Long Tab & 3" x 1.25" Base

**HCPCS:** A4333 / A5200

The CATH•SECURE®, from the M.C. Johnson Company, features an extended tab to secure larger sized tubes. It provides an efficient, standardized system for securing Foley catheters and various other medical tubing without restricting the flow. The device is easy to adjust and may remain in place for approximately one week without altering skin integrity. The open-close tab permits repeated access to the tubing, reducing medical tube trauma and eliminating the need for constant retaping.

## Key Features and Benefits

### CATH•SECURE® EXTENDED TAB

- **Low Profile;** less likely to get caught on clothing or snag.
- **Multiple Applications;** SKU consolidation, can be used for all types of tubing.
- **Easy to apply, easy to adjust;** velcro tab allows for easy access to tubing, eliminating the need for constant retaping.
- **Longevity;** device may stay in place for approximately one week without altering skin integrity.
- **Made in America**
- **DEHP FREE, HYPOALLERGENIC, NOT MADE WITH NATURAL RUBBER LATEX**

## Uses

Indwelling Urinary Catheters  
Central Line Ports  
Gastrostomy  
Jejunostomy  
Peritoneal Dialysis Tubing  
Hemovac Wound Drainage

Jackson-Pratt Wound Drainage  
EKG Lead Wires  
Prostatectomy Tubes  
Leg Bag Tubing  
Condom Catheters  
PICC Lines



M. C. Johnson Company, Inc.

8801 BUSINESS PARK DRIVE | SUITE 201 | FORT MYERS | FLORIDA | 33912  
239.433.0619 | 800.248.4058 | fax 239.433.1498 | www.mcjohnson.com