



Born of
Necessity

snappies

Collect,
store and
feed using
a single
innovative
breast milk
container





Leak-resistant seal

Flip-top hinge

Built in label section

Graduations in
both mL and oz.

Stable base



COLLECT

Aseptic and Ready to Use: Manufactured with an aseptic interior.

Minimize Contamination: The flip-top is conveniently hinged to the container so it will never be lost or dropped and can be easily opened and closed with one hand.

Convenient: The lid locks back out of the way while in use and closes with an audible “snap” to indicate an airtight, leak-resistant seal.

Direct connection: Compatible with Ameda®, Medela® and most other standard ring-sized pumps.



STORE

Leak Proof: Our patented triple-seal design provides an airtight and leak-resistant seal for maximum shelf life in a refrigerator or freezer

Easy Identification: Record the mother’s name, date and collection time to eliminate any potential identification problems.



Ideal Material: Food-grade polypropylene will not bind or hold milk fat during storage. Does not contain Bisphenol A, PVC or DEHP.

Stable: Large base keeps Snappies stable container while in use or when stored on a shelf in refrigerator or freezer.



FEED

No Transfer: The baby can be fed from the same container by mothers and other caregivers.

Convenience Sizes: 1oz. (35mL) is ideal for colostrum collection. 2.3 oz. (70mL) size is ideal for NICU and young babies.

Dual Graduations: Accurate, easy-to-read graduations in both oz. and mL.

Disposable and Recyclable: Making the world better for the next generation.



1 oz. Snappies

- Ideal for colostrum collection – minimize loss of ‘Precious Gold’
- Can attach directly to breastpump or can hand express into
- Appropriately sized to store mother’s milk as her supply increases to meet baby’s needs
- Requires less space for storage
- Ideal size for feeding, minimizing waste
- Nursing staff can easily recognize the smaller containers and use them first

2.3 oz. Snappies

- Ideal when Mother’s milk supply is fuller
- Ideal quantity for NICU baby feeds – no waste of precious breastmilk

11 mL Colostrum Collector

- Colostrum collectors minimize loss of ‘Precious Gold’
- Ideal for collection of hand expressed colostrum
- Sometimes it is good to hand express colostrum prior to delivery
- After collection store in the freezer
- When it is time to deliver, bring the collected colostrum to the hospital or birth center

The aseptic, all-in-one solution to simplify the collection, storage and feeding of expressed breast milk. Eliminates the need to transfer precious breast milk between containers and thus lowers the chance of contamination and risk of infection.

For further information and to order Snappies containers, please contact 800-772-8871 visit us online at www.snappiescontainers.com

**Snappies are BPA free,
disposable and recyclable.**



All Snappies products are proudly made in the USA.



snappies

**A leak-resistant, airtight solution
for the collection, storage and
feeding of breast milk.**

**Available in the 11 mL Colostrum
Collector, 1oz. and 2.3 oz. sizes!**

Used in hospitals and NIC units across North America, Snappies® enable moms to collect, store and feed breast milk from the same container. Our patented aseptic, airtight, leak-resistant containers give lactation consultants, care givers and mothers a helping hand to make breastfeeding easier in those early, critical months.



Ordering Information

Product	Packaging	Part #
1 oz. (35 mL)	10 per Heat-Sealed Bag, (100 total per case)	10256
2.3 oz. (70 mL)	10 per Heat-Sealed Bag, (200 total per case)	10253
2.3 oz. (70 mL)	Retail Pack, 12 per Box	10255
2.3 oz. (70 mL)	Retail Case, 12 Packs (144 total per case)	10259
11 mL Colostrum Collectors	Case (100 per case)	10275



2039 McMillan Street • Auburn AL 36832
Part of Thermo Fisher Scientific Inc.

For further information and to order Snappies containers,
please contact 800-772-8871
or visit us online at www.snappiescontainers.com

©2013 Thermo Fisher Scientific Inc. All rights reserved. Ameda is a registered trademark of Ameda AG. Medela is a registered trademark of Medela Inc. All other trademarks are the property of Thermo Fisher Scientific Inc. and its subsidiaries. Specifications, terms and pricing are subject to change. Not all products are available in all countries. Please consult your local sales representative for details.