



Alcohol Prep Pads

The Regard Alcohol Prep Pad is an absorbent, 2-ply pad, saturated in 70% isopropyl alcohol.

Our Alcohol Prep Pads:

- » Use foil packaging to maintain sterility and moisture
- » Are individually packed for single-use
- » Are made from a 2-ply, non-woven material



Products Available

REORDER #	DESCRIPTION	SIZE	PUOM	PUOM FACTOR	PACK STRING
3150831	Alcohol Prep Pad, 2-Ply	M	CA	4,000	200 packs/box, 20 boxes/ca

General Information

Active Ingredient: Isopropyl Alcohol (70%)

Pad Material: Laminated Non-Woven

Ply: 2-Ply

Pad Size (folded): 30mm x 30mm

Pad Size (unfolded): 30mm x 60mm

Pouch Size: 50mm x 50mm

Package: Aluminum Foil Film

Single-use: Yes

Sterile: Yes

Made w/ Natural Rubber Latex: No

Color: White

Standard: %IPA, Liquid Gram Weight

Product Facts

Use: Aids in the removal of oils and residue.

Warnings: Keep out of reach of children. In case of accidental ingestion, seek professional assistance or consult a poison control center immediately.

Directions: Apply locally as needed.

Ingredients: Isopropyl Alcohol 70%, purified water.

ORDERING & SERVICE INFO

support@roiscs.com

1-855-845-9302

Specifications subject to change without prior notice.

IMPORTANT ORDERING INSTRUCTIONS

Orders must be placed directly with ROi using ROi reorder numbers to receive HealthTrust contracted pricing.

Specifications subject to change without prior notice.



Distributed by
Resource Optimization
& Innovation, LLC
Nashville, TN 37203