



Anti-Fog Surgical Mask, Tie

ASTM Level 3

The Regard ASTM Level 3 Anti-Fog Surgical Mask provides maximum barrier protection in environments with a high risk of fluid exposure. The Anti-Fog technology prevents condensation from escaping the top of the mask, leaving glasses clear.



Our Anti-Fog Surgical Masks are:

- » Made with ties designed to provide universal size
- » Embossed with the protection level on the outside for easy identification
- » Suitable for a broad range of surgical applications, including laser procedures (refer to AORN and applicable standards for compliance)



Anti-Fog Surgical Mask, Tie

ASTM Level 3

ASTM Level 3 provides maximum barrier protection and is suitable for use in environments with a high risk of fluid exposure.

Products

REORDER #	DESCRIPTION	PUOM	PUOM FACTOR	PACK STRING
7757074	Anti-Fog Surgical Mask w/ Tie, ASTM L3	CA	1,000	50 ea/box, 20 boxes/ca

General Information

BFE/PFE: 99.5%

Material: Polypropylene (PP) Spunbond,
Polypropylene (PP) Meltblown

Style: 3-Ply, pleated with Tie / Anti-Fog Strip

Sterile: No

Made w/ Natural Rubber Latex: No

Color: Blue

Standard: ASTM F2100

ORDERING & SERVICE INFO

support@roiscs.com
1-855-845-9302

Specifications subject to change without prior notice.

IMPORTANT ORDERING INSTRUCTIONS

Orders must be placed directly with ROi using ROi reorder numbers to receive HealthTrust contracted pricing. Specifications subject to change without prior notice.



Distributed by
Resource Optimization
& Innovation, LLC
Nashville, TN 37203