



# Behavioral Health Mask, Earloop

## ASTM Level 3

The Regard ASTM Level 3 Behavioral Health Mask is an excellent choice for behavioral health patients, patient exams and visitors, and healthcare staff. The soft, lightweight material is breathable, yet provides maximum barrier protection against airborne particles and is suitable for use in environments with a high risk of fluid exposure.

---

### Our Behavioral Health Masks are:

- » Designed with elastic earloops
- » Universally sized
- » Made without a nose adjuster for behavioral health



No nose wire for use in a setting where metal is unacceptable.



# Behavioral Health Mask, Earloop

## ASTM Level 3

*ASTM Level 3 provides maximum barrier protection and is suitable for use in environments with a high risk of fluid exposure.*

### Products

REORDER #	DESCRIPTION	PUOM	PUOM FACTOR	PACK STRING
7700908	Behavioral Mask, ASTM L3, Earloop	CA	1,000	50 ea/box, 20 boxes/ca

### General Information

**BFE/PFE:** 99.5%

**Material:** Polypropylene (PP) Spunbond,  
Polypropylene (PP) Meltblown

**Style:** 3-Ply, Pleated with Earloops

**Sterile:** No

**Made w/ Natural Rubber Latex:** No

**Color:** Green

**Standard:** ASTM F2100

#### ORDERING & SERVICE INFO

support@roiscs.com  
1-855-845-9302

Specifications subject to change without prior notice.

#### IMPORTANT ORDERING INSTRUCTIONS

Orders must be placed directly with ROi using ROi reorder numbers to receive HealthTrust contracted pricing. Specifications subject to change without prior notice.



Distributed by  
Resource Optimization  
& Innovation, LLC  
Nashville, TN 37203