



Bouffant Cap

The Regard Bouffant Caps are breathable and lightweight while still providing a comfortable and secure fit.

Our Bouffant Caps are:

- » Available in two material types:
non-woven Spunbond
Polypropylene and Spunbond
Meltblown Spunbond
- » Elastic to provide a comfortable and secure fit without binding
- » Available in multiple colors





Bouffant Cap

Products

REORDER #	DESCRIPTION	SIZE	COLOR	PUOM	PUOM FACTOR	PACK STRING
7545092	Bouffant Cap, Polypropylene	24in	Red	CA	500	100 ea/box, 5 boxes/ca
7545087	Bouffant Cap, SMS	24in	Blue	CA	500	100 ea/box, 5 boxes/ca
3137387	Bouffant Cap, Polypropylene	24in	Blue	CA	500	100 ea/box, 5 boxes/ca
7545103	Bouffant Cap, SMS	28in	Blue	CA	500	100 ea/box, 5 boxes/ca

General Information

Material: SPP, SMS

Elastic: Polyamide

Seams: Sewn

Size: Universal

Fluid Resistance: No

Disposable: Yes


Sterile: No

Made w/ Natural Rubber Latex: No

Colors: Blue, Red

Sizing Guide (unit: in)

REORDER #	SIZE	A (RELAXED)	B (STRETCHED)
7545092	24"	7.1	22.5
7545087	24"	7.1	22.5
7545103	28"	10.6	27.0
3137387	24"	8.7	24.3
	Tolerance	±0.8	±0.8



ORDERING & SERVICE INFO

support@roiscs.com
1-855-845-9302

Specifications subject to change without prior notice.

IMPORTANT ORDERING INSTRUCTIONS

Orders must be placed directly with ROi using ROi reorder numbers to receive HealthTrust contracted pricing. Specifications subject to change without prior notice.



Distributed by
Resource Optimization
& Innovation, LLC
Nashville, TN 37203