



Disposable Scrub Shirts

The Regard Disposable Scrub Shirts provide a comfortable apparel choice for patients and staff. The V-neck design of the shirt makes donning and doffing comfortable and easy.

Our Disposable Scrub Shirts are:

- » Made of soft, anti-static SMMS for an added layer of protection
- » Designed with lightweight, breathable materials for added comfort
- » Constructed with three pockets for convenience





Disposable Scrub Shirts

Products

REORDER #	DESCRIPTION	SIZE	PUOM	PUOM FACTOR	PACK STRING
7567298	Scrub Shirt, Disposable	M	CA	94	94/ca (47 pcs/bag, 2 bags/ca)
7567302	Scrub Shirt, Disposable	L	CA	84	84/ca (42 pcs/bag, 2 bags/ca)
7567301	Scrub Shirt, Disposable	XL	CA	74	74/ca (37 pcs/bag, 2 bags/ca)
7567297	Scrub Shirt, Disposable	2XL	CA	74	74/ca (37 pcs/bag, 2 bags/ca)

General Information

Material: Soft SMMS Anti-static

Neck Style: V-Neck

Pockets: Yes

Sleeve Seal: Sewn (3 threads)

Water Vapor Transmission Rate: >3700 (single)

Disposable: Yes

Sterile: No


Made w/ Natural Rubber Latex: No

Color: Blue

Standard: ISO9073-3; DIN53923

Sizing Guide (unit: cm)

REF		M	L	XL	2XL	Tolerance
A	Shirt length	74.0	76.0	78.0	80.0	± 1.5
B	Side length (bottom to arm curvature)	34.0	35.0	36.0	38.0	± 2.0
C	Arm opening width	20.0	21.0	21.5	22.5	± 1.0
D	Neck width	20.3	23.3	23.3	24.0	± 1.5
E	V-neck length	14.6	14.6	14.6	14.6	± 1.0
F	Shirt bottom width	57.0	61.0	65.0	70.0	± 1.5
G	Front pocket width	17.0	17.0	17.0	17.0	± 0.5
H	Front pocket length	19.5	19.5	19.5	19.5	± 0.5
I	Breast pocket width	15.0	15.0	15.0	15.0	± 0.5
J	Breast pocket length	17.5	17.5	17.5	17.5	± 0.5
K	Shoulder seam length	30.5	31.0	33.0	36.5	± 1.0



ORDERING & SERVICE INFO

support@roiscs.com
1-855-845-9302

Specifications subject to change without prior notice.

IMPORTANT ORDERING INSTRUCTIONS

Orders must be placed directly with ROi using ROi reorder numbers to receive HealthTrust contracted pricing. Specifications subject to change without prior notice.



Distributed by
Resource Optimization
& Innovation, LLC
Nashville, TN 37203



Disposable Scrub Pants

The Regard Disposable Scrub Pants provide a comfortable apparel choice for patients and staff. The elastic waist design provides a secure and comfortable fit.

Our Disposable Scrub Pants are:

- » Made of soft, anti-static SMMS for an added layer of protection
- » Designed with lightweight, breathable materials for added comfort
- » Constructed with two front pockets for convenience





Disposable Scrub Pants

Products

REORDER #	DESCRIPTION	SIZE	PUOM	PUOM FACTOR	PACK STRING
7310425	Scrub Pants, Disposable	M	CA	70	70/ca (35 pcs/bag, 2 bags/ca)
7567300	Scrub Pants, Disposable	L	CA	70	70/ca (35 pcs/bag, 2 bags/ca)
7567296	Scrub Pants, Disposable	XL	CA	60	60/ca (30 pcs/bag, 2 bags/ca)
7567299	Scrub Pants, Disposable	2XL	CA	60	60/ca (30 pcs/bag, 2 bags/ca)

General Information

Material: Soft SMMS Anti-static

Pockets: Yes

Waist: Elastic

Pant Seal: Sewn (3 threads)

Water Vapor Transmission Rate: >3700 (single)

Disposable: Yes

Sterile: No


Made w/ Natural Rubber Latex: No

Color: Blue

Standard: ISO9073-3; DIN53923

Sizing Guide (unit: cm)

REF		M	L	XL	2XL	Tolerance
A	Hip opening, relaxed	27.5	31.0	33.0	35.0	± 3.0
A	Hip opening, fully stretched	55.0	59.0	63.0	67.0	± 3.0
B	Pant length	106.0	108.0	110.0	116.0	± 3.0
C	Crotch length	30.0	31.0	31.5	32.5	± 1.5
D	Leg bottom opening width	23.5	24.0	24.5	25.0	± 2.0
E	Pocket patch exterior width	17.0	17.0	17.0	17.0	± 1.0
F	Pocket length, leg side	16.0	16.0	16.0	16.0	± 1.0
G	Pocket length, crotch side	28.0	28.0	28.0	28.0	± 1.0
H	Upper thigh width	34.5	37.0	39.5	42.0	± 2.0
I	Inseam length	74.5	76.5	78.5	83.5	± 2.0



ORDERING & SERVICE INFO

support@roiscs.com
1-855-845-9302

Specifications subject to change without prior notice.

IMPORTANT ORDERING INSTRUCTIONS

Orders must be placed directly with ROi using ROi reorder numbers to receive HealthTrust contracted pricing. Specifications subject to change without prior notice.



Distributed by
Resource Optimization
& Innovation, LLC
Nashville, TN 37203