



Emesis Bags & Bag Dispensers

The Regard Emesis Bags and Bag Dispenser feature a twist and hook system that prevents spills and controls odors, making it easy to dispose of contents safely.



Our Emesis Bags & Bag Dispensers feature:

- » A rigid Polypropylene ring that provides a wide opening to ensure unobstructed access to the containment bag
- » A cardboard rim option with slim packaging, designed for convenient storage
- » Emesis Bags that are constructed of tear-resistant HDPE or LDPE and can safely hold up to 36 fluid ounces
- » An easy-to-load dispenser that can securely hold 24 rigid ring emesis bags
- » Anchors and screws for convenient wall mounting





Emesis Bags

Products

REORDER #	DESCRIPTION	SIZE	BAG MATERIAL	COLOR	PUOM	PUOM FACTOR	PACK STRING
7611844	Emesis Bag w/ Rigid Ring & Graduations	36oz/ Max 1500ml	HDPE	Green	CA	144	24 ea/bag, 6 bags/ca
7611845	Emesis Bag w/ Rigid Ring & Graduations	36oz/ Max 1500ml	LDPE	Blue	CA	144	24 ea/bag, 6 bags/ca
9008851	Emesis Bag w/ Cardboard Rim	36oz	LDPE	Clear	CA	48	1 ea/bag, 48 bags/ca

General Information

Bag Material: HDPE, LDPE

Ring Material: Polypropylene, Paperboard (9008851)

Bag Size: 36oz, max fill 1500ml (7611844, 7611845); 36oz (9008851)

Graduated Markings: Every 4 fl oz (100cc/100ml) up to 36 fl oz (1000cc/1000ml) with 7611844, 7611845 max fill 1500cc/1500ml

Colors: Blue (7611844), Green (7611845), Clear (9008851)

Single-use: Yes

Sterile: No

Made w/ Natural Rubber Latex: No

Compatible Dispenser: 7611846 for 7611844, 7611845

Emesis Bag Dispenser

Products

REORDER #	DESCRIPTION	PUOM	PUOM FACTOR	PUOM STRING
7611846	Dispenser Box, Round	CA	6	1/box, 6 boxes/ca

General Information

Material: Acrylic, Silicon

Shape: Round

Capacity: 24 Regard Emesis Bags

Body Color: Transparent

Lid Color: Blue

Single-use: No

Sterile: No

Made w/ Natural Rubber Latex: No

Compatible Emesis Bags: 7611844, 7611845

ORDERING & SERVICE INFO

support@roiscs.com
1-855-845-9302

Specifications subject to change without prior notice.

IMPORTANT ORDERING INSTRUCTIONS

Orders must be placed directly with ROi using ROi reorder numbers to receive HealthTrust contracted pricing. Specifications subject to change without prior notice.



Distributed by
Resource Optimization
& Innovation, LLC
Nashville, TN 37203