



Heavy-Duty Knee-High Boot Covers

with Ankle Ties

The Regard Heavy-Duty Knee-High Boot Covers are designed with Polypropylene Laminate material to offer optimal barrier protection from knee to toe.

Our Heavy-Duty Knee-High Boot Covers:

- » Are made from fluid-resistant laminate for durability, breathability and comfort
- » Feature an elastic band at the top and ankle ties for fit, support, and protection
- » Are designed with non-skid treads for added safety in slippery conditions





Heavy-Duty Knee-High Boot Covers with Ankle Ties

Products

REORDER #	DESCRIPTION	SIZE	PUOM	PUOM FACTOR	PACK STRING
7545089	Knee-High Boot Cover w/ Elastic Band Top & Poly-coated Tie Closure	L	CA	120	30 ea/box, 4 boxes/ca
7545098	Knee-High Boot Cover w/ Elastic Band Top & Poly-coated Tie Closure	XL	CA	120	30 ea/box, 4 boxes/ca

General Information

Material: Laminate

Non-skid: Foam

Seams: Ultrasonic

Fluid Resistance: Yes

Disposable: Yes


Sterile: No

Made w/ Natural Rubber Latex: No

Color: Blue

Sizing Guide (unit: cm)

REORDER #	SIZE	A	B	C	D
7545089	L	37.0	28.5	29.0	30.0
7545098	XL	40.0	28.5	31.5	35.5
	Tolerance	±1.5	±0.75	±1.5	±1.5



ORDERING & SERVICE INFO

support@roiscs.com

1-855-845-9302

Specifications subject to change without prior notice.

IMPORTANT ORDERING INSTRUCTIONS

Orders must be placed directly with ROi using ROi reorder numbers to receive HealthTrust contracted pricing. Specifications subject to change without prior notice.



Distributed by
Resource Optimization
& Innovation, LLC
Nashville, TN 37203