



Heavy-Duty Shoe Covers

The Regard Heavy-Duty Shoe Cover is made with extra strong fluid resistant SMS material for increased durability.

Our Heavy-Duty Shoe Covers are:

- » Made from Low-Lint, Spunbond Meltblown Spunbond for breathability, strength and comfort
- » Designed with non-skid treads for safety in slippery conditions
- » Constructed with ultrasonic welded inseams





Heavy-Duty Shoe Covers

Products


REORDER #	DESCRIPTION	SIZE	PUOM	PUOM FACTOR	PACK STRING
7545088	Ankle-High Shoe Cover w/ Elastic Band	LG	CA	300	100 ea/box, 3 boxes/ca
7545085	Ankle-High Shoe Cover w/ Elastic Band	XL	CA	240	80 ea/box, 3 boxes/ca
7545095	Ankle-High Shoe Cover w/ Elastic Band	2XL	CA	400	50 ea/box, 8 boxes/ca

General Information

Material: SMS
Non-skid: Yes
Seams: Ultrasonic Welded
Fluid Resistance: Yes
Disposable: Yes
Sterile: No
Made w/ Natural Rubber Latex: No
Color: Blue

Sizing Guide (unit: cm)

REORDER #	SIZE	A	B	C	D	E
7545088	LG	37.0	19.0	28.0	26.0	1.5
7545085	XL	41.0	21.5	31.0	26.5	1.5
7545095	2XL	44.5	21.5	34.5	29.5	1.5
	Tolerance	±1.3	±1.3	±1.3	±1.3	±0.5



ORDERING & SERVICE INFO

support@roiscs.com
1-855-845-9302

Specifications subject to change without prior notice.

IMPORTANT ORDERING INSTRUCTIONS

Orders must be placed directly with ROi using ROi reorder numbers to receive HealthTrust contracted pricing. Specifications subject to change without prior notice.



Distributed by
Resource Optimization & Innovation, LLC
 Nashville, TN 37203