



# Air Assisted Lateral Transfer Mattresses, Kits & Accessories

## in collaboration with HoverTech

Designed to be multi-functional, the breathable Regard Air Assisted Lateral Transfer Mattress, powered by HoverTech, makes patient transfers, boosting, turning, positioning and proning easier, while protecting caregiver safety.



### Our Air Assisted Lateral Transfer Mattresses:

- » Comfortably cradle the patient in a stable position for safer patient handling
- » Enable caregiver to safely transfer patient without lifting or straining
- » Reduce the force required to move a patient by 80-90%





# Air Assisted Lateral Transfer Mattresses, Kits & Accessories

## Products

REORDER #	DESCRIPTION	SIZE	WEIGHT LIMIT	PUOM	PUOM FACTOR	PACK STRING
7629489	Air Assisted Full-Size Matt	34" W x 78" L	1,200 lbs	CA	10	10/ca
7629490	Air Assisted Full-Size Matt	39" W x 78" L	1,200 lbs	CA	10	10/ca
7629491	Air Assisted Full-Size Matt	50" W x 78" L	1,200 lbs	CA	5	5/ca
7629492	Air Assisted Half-Size Matt	34" W x 45" L	600 lbs	CA	10	10/ca
7629493	Air Assisted Half-Size Matt	39" W x 45" L	600 lbs	CA	10	10/ca

## Kits & Accessories

REORDER #	DESCRIPTION	DETAILS	PUOM	PUOM FACTOR	PACK STRING
8997499	Air Repositioning & Offloading Kit	39" Regard Matt, Full HoverCover, PROSWedge Pair with & without Tail	CA	5	5/ca
9004780	Air Repositioning & Offloading Kit	39" Regard Matt and Full HoverCover	CA	10	10/ca
2167648	Half HoverCover, 40" W x 40" L	Disposable Breathable Absorbent Cover	CA	50	50/ca
1743996	Full HoverCover, 39" W x 65" L	Disposable Breathable Absorbent Cover	CA	50	50/ca
7861756	PROSWedges Pair	One with Tail and One without Tail	CA	5	5/ca

## Blended Codes Program

The Blended Code Program offers blended codes that consist of new and reprocessed products. Customers need to participate in the Cardinal Health Used Product Collection Program to access the blended codes.

For more information contact Customer Support at [855-845-9302](tel:855-845-9302) or Email [support@roiscs.com](mailto:support@roiscs.com).

REORDER #	DESCRIPTION	SIZE	WEIGHT LIMIT	PUOM	PUOM FACTOR	PACK STRING
8981027	Air Assisted Full Size Matt	34" W x 78" L	1,200 lbs	CA	10	10/ca
8981026	Air Assisted Full Size Matt	39" W x 78" L	1,200 lbs	CA	10	10/ca
8981025	Air Assisted Full Size Matt	50" W x 78" L	1,200 lbs	CA	5	5/ca
8981024	Air Assisted Half Size Matt	34" W x 45" L	600 lbs	CA	10	10/ca
8981023	Air Assisted Half Size Matt	39" W x 45" L	600 lbs	CA	10	10/ca

## General Information

**Material:** Non-Woven Polypropylene Fiber (top), Nylon Fabric (bottom)

**Reprocessed Units Decontaminated:** Yes; w/ Ethylene Oxide

**Seams:** Sewn

**Sterile:** No

**Color:** White, Blue

**Made w/ Natural Rubber Latex:** No

## Equipment Procurement Options

The HoverTech Air Supply Equipment is compatible with all Regard Air Transfer Mattresses. It can be purchased direct from HoverTech for a separate fee or placed with customer through ROi's Placed Equipment Program.

The Placed Equipment Program requires execution of a separate agreement with ROi and disposable commitment. For more information on this program, please contact customer service at [support@roiscs.com](mailto:support@roiscs.com) or call [855-845-9302](tel:855-845-9302).