



Procedure Mask, Earloop

ASTM Level 1 & Level 3

The Regard Disposable ASTM Level 1 and Level 3 Procedure Masks with earloop design are flexible to universally fit staff with the protection level embossed on the outside of the mask for easy identification.

Our Procedure Masks are:

- » Designed with elastic earloops
- » Universally sized
- » Soft and comfortable to wear





Procedure Mask, Earloop

ASTM Level 1 & Level 3

ASTM Level 1 is recommended for use as a fluid-barrier solely in low-risk environments. ASTM Level 3 provides maximum barrier protection and is suitable for use in environments with a high risk of fluid exposure.

Products

REORDER #	DESCRIPTION	PUOM	PUOM FACTOR	PACK STRING
3137388	Procedure Mask w/ Earloop, ASTM L1	CA	1,000	50 ea/box, 20 boxes/ca
3137773	Procedure Mask w/ Earloop, ASTM L3	CA	1,000	50 ea/box, 20 boxes/ca

General Information

BFE/PFE: 99.5%

Material: Polypropylene (PP) Spunbond,
Polypropylene (PP) Meltblown

Style: 3-Ply, Pleated with Earloops

Sterile: No

Made w/ Natural Rubber Latex: No

Color: Green

Standard: ASTM F2100

ORDERING & SERVICE INFO

support@roiscs.com
1-855-845-9302

Specifications subject to change without prior notice.

IMPORTANT ORDERING INSTRUCTIONS

Orders must be placed directly with ROi using ROi reorder numbers to receive HealthTrust contracted pricing. Specifications subject to change without prior notice.



Distributed by
Resource Optimization
& Innovation, LLC
Nashville, TN 37203