



Standard Knee-High Boot Covers

with Elastic Bands

The Regard Standard Knee-High Boot Covers are designed to provide protection from knee to toe.

Our Standard Knee-High Boot Covers:

- » Are made from fluid-resistant Laminate and Low-Lint, Spunbond Meltblown Spunbond for durability, breathability and comfort
- » Feature an elastic band at the top and ankle for fit, support and protection
- » Are designed with non-skid treads for added safety in slippery conditions





Standard Knee-High Boot Covers with Elastic Bands

Products

REORDER #	DESCRIPTION	SIZE	PUOM	PUOM FACTOR	PACK STRING
7545094	Knee-High Boot Cover w/ Elastic Band Top & Ankle	L	CA	150	50 ea/box, 3 boxes/ca
7545086	Knee-High Boot Cover w/ Elastic Band Top & Ankle	XL	CA	150	50 ea/box, 3 boxes/ca


General Information

Material: Laminate, SMS
Non-skid: Tread
Seams: Heat-sealing
Fluid Resistance: Yes

Disposable: Yes
Sterile: No
Made w/ Natural Rubber Latex: No
Color: Blue

Sizing Guide (unit: cm)

REORDER #	SIZE	A	B	C	D	E
7545094	L	35.5	23.0	28.0	37.0	28.0
7545086	XL	42.5	23.0	32.0	38.5	29.0
	Tolerance	±1.5	±1.5	±1.5	±1.5	±1.5



ORDERING & SERVICE INFO

support@roiscs.com
1-855-845-9302

Specifications subject to change without prior notice.

IMPORTANT ORDERING INSTRUCTIONS

Orders must be placed directly with ROi using ROi reorder numbers to receive HealthTrust contracted pricing. Specifications subject to change without prior notice.



Distributed by
Resource Optimization
& Innovation, LLC
Nashville, TN 37203