



Standard Shoe Covers

The Regard Standard Shoe Cover is designed to be durable while maintaining a comfortable universal fit.

Our Standard Shoe Covers are:

- » Made from Low-Lint, Spunbond Polypropylene for breathability, strength and comfort
- » Designed with non-skid treads for safety in slippery conditions
- » Constructed with stitched inseams





Standard Shoe Covers

Products


REORDER #	DESCRIPTION	SIZE	PUOM	PUOM FACTOR	PACK STRING
7545099	Ankle-High Shoe Cover	LG	CA	400	100 ea/box, 4 boxes/ca
7545101	Ankle-High Shoe Cover	XL	CA	400	100 ea/box, 4 boxes/ca
7545097	Ankle-High Shoe Cover	2XL	CA	400	100 ea/box, 4 boxes/ca

General Information

Material: SPP **Disposable:** Yes
Non-skid: Yes **Sterile:** No
Seams: Sewn **Made w/ Natural Rubber Latex:** No
Fluid Resistance: No **Color:** Blue

Sizing Guide (unit: cm)

REORDER #	SIZE	A	B	C	D	E	F
7545099	LG	22.3	17.8	36.2	15.2	5.1	1
7545101	XL	23.5	19.7	40.0	15.2	5.1	1
7545097	2XL	24.8	20.3	41.3	15.9	5.1	1
	TOLERANCE	+2.5/-0	+2.5/-0.6	+1.2/-0	+2.5/-0.6	+1.2/-0.5	+0/-1



ORDERING & SERVICE INFO

support@roiscs.com
1-855-845-9302

Specifications subject to change without prior notice.

IMPORTANT ORDERING INSTRUCTIONS

Orders must be placed directly with ROi using ROi reorder numbers to receive HealthTrust contracted pricing. Specifications subject to change without prior notice.



Distributed by
Resource Optimization
& Innovation, LLC
Nashville, TN 37203