



# Surgeon Cap

## with Elastic Back

The Regard Surgeon Cap with Elastic Back offers protection against particulate while maintaining a secure fit with the elastic back.

---

### Our Surgeon Caps with Elastic Backs are:

- » Made from Spunbond Polypropylene and Spunbond Meltblown Spunbond material
- » Lightweight with a comfortable fit
- » Easy to don





# Surgeon Cap with Elastic Back

## Products

REORDER #	DESCRIPTION	SIZE	PUOM	PUOM FACTOR	PACK STRING
7545091	Surgeon Cap w/ Elastic Back	Universal	CA	500	100 ea/box, 5 boxes/ca

## General Information

**Material:** SPP, SMS

**Seams:** Sewn

**Back Fixation:** Elastic

**Size:** Universal

**Fluid Resistance:** No

**Disposable:** Yes


**Sterile:** No

**Made w/ Natural Rubber Latex:** No

**Color:** Blue

## Sizing Guide (unit: cm)

REORDER #	SIZE	A	B	C	D (STRETCHED)	E	F (ELASTIC RELAX)
7545091	Universal	22.9	17.1	13.6	31.9	19.1	10.2
	<b>TOLERANCE</b>	±1.3	±1.3	±1.3	±2	±1.3	±0.5



### ORDERING & SERVICE INFO

support@roiscs.com

1-855-845-9302

Specifications subject to change without prior notice.

### IMPORTANT ORDERING INSTRUCTIONS

Orders must be placed directly with ROi using ROi reorder numbers to receive HealthTrust contracted pricing.

Specifications subject to change without prior notice.



Distributed by  
Resource Optimization  
& Innovation, LLC  
Nashville, TN 37203