



Surgeon Hood with Beard Cover

The Regard Surgeon Hood with Beard Cover provides cool and comfortable coverage of the head and neck.

Our Surgeon Hoods with Beard Covers are:

- » Made from Spunbond Polypropylene and Spunbond Meltblown Spunbond
- » Adjustable with neck ties for maximum coverage of hair and sideburns
- » Universally sized





Surgeon Hood with Beard Cover

Products


REORDER #	DESCRIPTION	SIZE	COLOR	PUOM	PUOM FACTOR	PACK STRING
7545093	Surgeon Hood & Beard Cover w/ Under Chin Tie	Universal	Blue	CA	300	100 ea/box, 3 boxes/ca

General Information

Material: SPP, SMS **Disposable:** Yes
Beard Cover: Yes **Sterile:** No
Seams: Sewn **Made w/ Natural Rubber Latex:** No
Size: Universal **Color:** Blue
Fluid Resistance: No

Sizing Guide (unit: cm)

REORDER #	SIZE	A	B	C	D	E	F	G	H	I
7545093	Universal	26.7	16.8	5.1	25.2	31.1	23.5	71.1	8.3	45.7
	TOLERANCE	±1	±1	±0.5	±1.5	±1.5	±1.5	±1.5	±0.5	±1.5



ORDERING & SERVICE INFO

support@roiscs.com
1-855-845-9302

Specifications subject to change without prior notice.

IMPORTANT ORDERING INSTRUCTIONS

Orders must be placed directly with ROi using ROi reorder numbers to receive HealthTrust contracted pricing. Specifications subject to change without prior notice.



Distributed by
Resource Optimization & Innovation, LLC
 Nashville, TN 37203