



Surgical Mask, Tie

ASTM Level 3

The Regard Surgical Mask is designed for a universal secure fit while still maintaining maximum barrier protection.

Our Surgical Masks are:

- » Made with ties designed to provide universal size
- » Easily adjusted for comfort and maximum protection
- » Breathable and fluid resistant





Surgical Mask, Tie

ASTM Level 3

ASTM Level 3 provides maximum barrier protection and is suitable for use in environments with a high risk of fluid exposure.

Products

REORDER #	DESCRIPTION	PUOM	PUOM FACTOR	PACK STRING
7307158	Surgical Mask w/Tie, ASTM L3	CA	1,000	50 ea/box, 20 boxes/ca

General Information

BFE/PFE: 99.5%

Material: Polypropylene (PP) Spunbond,
Polypropylene (PP) Meltblown

Style: 3-Ply, Pleated with Tie

Sterile: No

Made w/ Natural Rubber Latex: No

Color: Blue

Standard: ASTM F2100

ORDERING & SERVICE INFO

support@roiscs.com
1-855-845-9302

Specifications subject to change without prior notice.

IMPORTANT ORDERING INSTRUCTIONS

Orders must be placed directly with ROi using ROi reorder numbers to receive HealthTrust contracted pricing. Specifications subject to change without prior notice.



Distributed by
Resource Optimization
& Innovation, LLC
Nashville, TN 37203