



# Carter-Thomason Port Closure Systems


in collaboration with Cooper Surgical

Regard's Port Closure Portfolio provides a comprehensive offering to ensure consistent suture passings on any body type. Designed for safe and efficient port site closure for both Laparoscopic and Robotic procedures including Ventral Hernia Repair.

---

## Our Carter-Thomason Port Closure Systems feature:

- » Ergonomically engineered handles
- » Rigid needle to enable precision control
- » The choice of two distinct Guide/Passer designs





# Carter-Thomason Port Closure Systems



Reorder # 733776

## Carter-Thomason CloseSure System XL

**Single-use kit includes:**

- (1) 10/12 mm Pilot Guide
- (1) 15 mm Pilot Guide
- (1) XL CloseSure Suture Passer



Reorder # 75926

## Carter-Thomason CloseSure System

**Single-use kit includes:**

- (1) 5 mm Pilot Guide
- (1) 10/12 mm Pilot Guide
- (1) CloseSure Suture Passer



Reorder # 775984

## Carter-Thomason II Port Closure System

**Single-use kit includes:**

- (1) 10/12 mm Suture Guide
- (1) 15 mm Suture Guide
- (1) SureGrasp Suture Passer



Reorder # 351277

## Carter-Thomason II Port Closure System

**Single-use kit includes:**

- (1) 10/12 mm Pilot Guide
- (1) 15 mm Pilot Guide
- (1) SureGrasp Suture Passer



Reorder # 534545

## Carter-Thomason II Suture Passer

**Single-use suture passer includes:**

- (1) SureGrasp Suture Passer

## Products

REORDER #	DESCRIPTION	PACK STRING
733776	Carter-Thomason CloseSure System XL, 10/12 mm & 15 mm Pilot Guides + XL Suture Passer	3 kits/box
75926	Carter-Thomason CloseSure System, 5 mm & 10/12 mm Pilot Guides + Suture Passer	5 kits/box
775984	Carter-Thomason II - Port Closure System, 10/12 mm & 15 mm Suture Guides + SureGrasp Suture Passer	5 kits/box
351277	Carter-Thomason II - Port Closure System, 10/12 mm & 15 mm Pilot Guides + SureGrasp Suture Passer	5 kits/box
534545	Carter-Thomason II - SureGrasp Suture Passer	5 each/box

## General Information

**Suture Passer Molded Body & Plunger:**

Polypropylene & Polycarbonate

**Suture Passer Shaft Material:** Stainless Steel

**Suture Passer Tip Geometry:** Beveled Needle Tip

**Suture Passer Diameter:** 0.028 cm

**Sterile:** Yes

**Color:** White



# IMPORTANT HealthTrust Members Ordering Instructions:

Regard has collaborated with Cooper Surgical to bring added value to their members through a co-branded Port Closure device. Under this co-branded model, HealthTrust members will now have exclusive access to the same best-in-class port closure devices they used before, with the added benefit of secured distribution and market-leading pricing.

Orders must be placed direct with ROi using ROi reorder numbers. Using any other number when ordering this product can lead to off-contract spend and higher costs.

ROI REORDER #	MANUFACTURER #	DESCRIPTION
733776	CTXL	Carter-Thomason CloseSure System XL, 10/12 mm & 15 mm Pilot Guides + XL Suture Passer
75926	CTI-512N	Carter-Thomason CloseSure System with 5 mm & 10/12 mm Pilot Guides + Suture Passer
775984	CTI-1015P	Carter-Thomason II - Port Closure System, 10/12 mm & 15 mm Suture Guides + SureGrasp Suture Passer
351277	CTI-1015SG	Carter-Thomason II - Port Closure System, 10/12 mm & 15 mm Pilot Guides + SureGrasp Suture Passer
534545	CTSG	Carter-Thomason II - SureGrasp Suture Passer

#### ORDERING & SERVICE INFO

support@roiscs.com  
1-855-845-9302

Specifications subject to change without prior notice.

#### IMPORTANT ORDERING INSTRUCTIONS

Orders must be placed directly with ROi using ROi reorder numbers to receive HealthTrust contracted pricing. Specifications subject to change without prior notice.



Distributed by  
Resource Optimization  
& Innovation, LLC  
Nashville, TN 37203



# Port Closure Label Information

Please note the only changes you may experience are labeling changes as highlighted below. The product itself is the same Carter-Thomason devices that your facilities have historically used.

## OFF-CONTRACT PRODUCT LABEL

Products with this **COOPER SURGICAL** label are off-contract and will result in higher pricing.

**PEEL HERE**

**CARTER-THOMASON**  
CloseSure System®

**REF** CTI-512N  
**LOT** XXXXXX  
YYYY-MM-DD

(01) 6 0888937 00009 4  
(10) XXXXXX  
(17) 000000  
(240) CTI-512N

**Contents:**  
(1) Carter-Thomason® Suture Passer  
(2) Pilot® Guides

**STERILE EO** **R<sub>x</sub> Only** **CE 2797**

- Do not reuse
- Do not resterilize
- Do not use if package is damaged
- Consult instructions for use
- Caution
- Not made with natural rubber latex

This product is protected under one or more of the following U.S. Patents:  
5,899,911 • 6,183,485

**CooperSurgical**  
Phone: (800) 243-2974 • Fax: (800) 262-0105  
International  
Phone: +1(203) 601-9818 • Fax: +1(203) 601-4747  
www.coopersurgical.com

**EC REP**  
EMERGO EUROPE  
Prinsessegracht 20  
2514 AP The Hague  
The Netherlands  
TL-CTI512N-C • Rev. A

CooperSurgical, Inc  
95 Corporate Drive  
Trumbull, CT 06611 USA  
Made in the USA

## ON-CONTRACT PRODUCT LABEL

Products with this **REGARD** label are on-contract and will result in best available pricing.

**PEEL HERE**

**CARTER-THOMASON**  
CloseSure System®

**REF** 75926  
**LOT** XXXXXX  
YYYY-MM-DD

(01) 6 0888937 02703 9  
(10) XXXXXX  
(17) 000000  
(240) 75926

**Contents:**  
One (1) Carter-Thomason® Suture Passer  
One (1) 5 mm Pilot® Guide  
One (1) 10/12 mm Pilot® Guide

**STERILE EO** **R<sub>x</sub> Only** **MD**

- Do not reuse
- Do not resterilize
- Do not use if package is damaged
- Consult instructions for use  
www.coopersurgical.com/en
- Caution
- Not made with natural rubber latex

**regard**

Resource Optimization & Innovation, LLC  
645 Maryville Center Drive  
St. Louis, Missouri 63141  
855-845-9302  
support@roiscs.com  
CLP4273461 • Rev. A

CooperSurgical®  
95 Corporate Drive  
Trumbull, CT 06611 USA  
Phone: (800) 243-2974  
Fax: (800) 262-0105  
Made in the USA