

DVT Therapy Portfolio





Kendall SCD™ 700

Smart Compression™ System

Regard® in collaboration with Cardinal Health offers Smart Compression, the next evolution of intermittent pneumatic compression (IPC). It's clinically designed and proven to move more blood¹ to help prevent stasis,¹ track compliance² and educate patients and clinicians about the risk of venous thromboembolism (VTE).



The five components of Smart Compression™ are:

- » Smart compression cycles
- » Smart sleeve design
- » Smart patient compliance
- » Smart sustainability
- » Smart training

1. Griffin M, Kakkos SK, Geroulakos G, Nicolaides AN. Comparison of three intermittent pneumatic compression systems in patients with varicose veins: a hemodynamic study. *Int Angiol.* 2007;26(2):158-164. 2. Tarone D. Selected long abstracts from the St. Luke's University Health Network Quality Awards Program. *Int J Acad Med.* 2017; 3: S176 -S18 8



Kendall SCD™ 700 Smart Compression™ System

Kendall SCD™ Comfort Sleeves

REORDER #	DESCRIPTION	SIZE	QUANTITY
393208	Knee Length	S	5 pairs/case
985992	Knee Length	M	5 pairs/case
985999	Knee Length	L	5 pairs/case
393214	Thigh Length	XS	5 pairs/case
575117	Thigh Length	S	5 pairs/case
985998	Thigh Length	M	5 pairs/case
985991	Thigh Length	L	3 pairs/case
986000	Thigh Length, Tear-Away	S	5 pairs/case
986001	Thigh Length, Tear-Away	M	5 pairs/case
986002	Thigh Length, Tear-Away	L	3 pairs/case



- Easy one-touch operation
- Compliance meter tracks patient's true wear time
- Eight-hour lithium ion battery
- Animated alerts show cause of error and how to resolve
- Durable design with IP23 rating
- Adjustable bed hook

Kendall SCD™ Express Comfort Sleeves

REORDER #	DESCRIPTION	SIZE	QUANTITY
375557	Knee Length	M	5 pairs/case
375963	Knee Length	L	5 pairs/case
166175	Knee Length	XL	5 pairs/case
375559	Thigh Length	S	5 pairs/case
375558	Thigh Length	M	5 pairs/case
375959	Thigh Length	L	3 pairs/case
409471	Thigh Length, Tear-Away	S	5 pairs/case
375550	Thigh Length, Tear-Away	M	5 pairs/case
411370	Thigh Length, Tear-Away	L	3 pairs/case
375940	Thigh Length, Sterile		5 pairs/case



- Circumferential, sequential and gradient
- Tear-Away Sleeve Technology optional
- Large and X-large sizes FDA cleared for use on bariatric patients
- Kendall SCD™ Comfort Sleeves incorporate Dry, Cool, Soft (DCS) Technology to address patient discomfort

Kendall SCD™ Express Foot Cuff

REORDER #	DESCRIPTION	SIZE	QUANTITY
375547	Foot Cuff	M	10 each/case
375546	Foot Cuff	L	10 each/case



- Optimize blood flow with Vascular Refill Detection
- Padded heel strap to minimize pressure points
- Fit left or right foot

Kendall SCD™ Tubing

REORDER #	DESCRIPTION	SIZE	QUANTITY
749797	Tubing Set	7 Feet	1 pair/box
75093	Tubing Set Extension	4 Feet	1 pair/box



- Tubing set provided with controller, seven feet
- Extension tubing set available, four feet



Kendall SCD™ 700

Smart Compression™ System

Blended Code Program

The Blended Code Program offers blended codes that consist of new and reprocessed products. Customers will need to participate in the Cardinal Used Product Collection Program to access the blended codes. For more information contact Customer Support at 855-845-9302 or Email Support@roiscs.com.



Consistent Value



620K pounds of medical waste diverted annually*

*Internal data available upon request.

Kendall SCD™ Comfort Sleeves

REORDER #	DESCRIPTION	SIZE	QUANTITY
2912116	Blended, Knee Length	S	25 pairs/case
2405397	Blended, Knee Length	M	25 pairs/case
2405398	Blended, Knee Length	L	25 pairs/case
754481	Blended, Thigh Length	XS	25 pairs/case
2405393	Blended, Thigh Length	S	25 pairs/case
7362937	Blended, Thigh Length	M	25 pairs/case
2405395	Blended, Thigh Length	L	15 pairs/case

Kendall SCD™ Express Comfort Sleeves

REORDER #	DESCRIPTION	SIZE	QUANTITY
1467943	Blended, Knee Length	M	25 pairs/case
2162378	Blended, Knee Length	L	25 pairs/case
2912121	Blended, Knee Length	XL	25 pairs/case
7362938	Blended, Thigh Length	S	25 pairs/case
2912118	Blended, Thigh Length	M	25 pairs/case
2912117	Blended, Thigh Length	L	15 pairs/case

Kendall SCD™ Express Foot Cuffs

REORDER #	DESCRIPTION	SIZE	QUANTITY
2912120	Blended, Foot Cuff	M	50 each/case
2912119	Blended, Foot Cuff	L	50 each/case

ORDERING & SERVICE INFO

support@roiscs.com
1-855-845-9302

Specifications subject to change without prior notice.

IMPORTANT ORDERING INSTRUCTIONS

Orders must be placed directly with ROi using ROi reorder numbers to receive HealthTrust contracted pricing. Specifications subject to change without prior notice.



Distributed by
Resource Optimization
& Innovation, LLC
Nashville, TN 37203



**SCAN THE QR CODE TO
DOWNLOAD THE REGARD CATALOG**

go.healthtrustpg.com/regard-catalog

www.roiscs.com • support@roiscs.com • 1.855.845.9302