



# Lubricating Jelly

Regard's Lubricating Jelly is formulated to spread evenly and increase patient comfort.

---

## Our Lubricating Jellies:

- » Are sterile and water soluble
- » Provide a smooth consistency
- » Are individually packed for single-use



# Lubricating Jelly

## Products Available

REORDER #	DESCRIPTION	SIZE	PUOM	PUOM FACTOR	PACK STRING
3150171	Lubricating Jelly	4oz Tube	CA	72	12 tubes/box, 6 boxes/ca
3149956	Lubricating Jelly	5gm Foil Pack	CA	1,500	150 packs/box, 10 boxes/ca

## General Information

**Active Ingredient:** Glycerin, Polyethylene Glycol-300, Polyethylene Glycol-1450, Carbomer, 4-Hydroxybenzoic Acid Methyl ester, 4-Hydroxybenzoic Acid Propylester, Sodium Hydroxide, Purified Water


**Package:** Aluminum Plastic Composite Tube (3150171)

**Storage Conditions:** Ambient Temperatures

**Single-use:** Yes

**Sterile:** Yes

**Made w/ Natural Rubber Latex:** No

<p><b>ORDERING &amp; SERVICE INFO</b> support@roiscs.com 1-855-845-9302 Specifications subject to change without prior notice.</p>	<p><b>IMPORTANT ORDERING INSTRUCTIONS</b> Orders must be placed directly with ROi using ROi reorder numbers to receive HealthTrust contracted pricing. Specifications subject to change without prior notice.</p>	 <p>Distributed by Resource Optimization &amp; Innovation, LLC Nashville, TN 37203</p>
--	---	---