



Surgical N95 Particulate Respirator

NIOSH Certified – ASTM Level 3 Fluid Resistant

Regard's Surgical N95 Particulate Respirator is a three-panel flat fold design used for at least 95 percent filtration efficiency against certain particulates. The unique design is comfortable to wear and fits a broad range of face sizes.



Our Surgical N95 Particulate Respirators are:

- » NIOSH certified
- » ASTM Level 3 Fluid Resistant
- » Fitted to form a seal around the nose and mouth
- » Equipped with a foam nose piece for enhanced comfort



Surgical N95 Particulate Respirator

NIOSH Certified – ASTM Level 3 Fluid Resistant

ASTM Level 3 provides maximum barrier protection and is suitable for use in environments with a high risk of fluid exposure.

Products

REORDER #	DESCRIPTION	PUOM	PUOM FACTOR	PACK STRING
7836620	Surgical N95 Particulate Respirator, NIOSH Certified, ASTM Level 3 Fluid Resistant	CA	400	20 ea/box, 20 boxes/ca

General Information

BFE/PFE: >98%

Material: Polypropylene (PP) Spunbond, Polyethylene Terephthalate (PET), Polypropylene (PP) Meltblown, Ethylene-Propylene Side By Side (ES), Polyester & Spandex (Strap)

Style: Flat Fold

Respirator Size: Standard

Fluid Resistance: Yes (160 mmHg)

Nose Wire Material: Aluminium

Sterile: No

Exhalation Valve: No

Made w/ Natural Rubber Latex: No

Cellulose Free: Yes

Color: White

Standard: ASTM F2100

NIOSH Approval Number: TC 84A-9745

ORDERING & SERVICE INFO

support@roiscs.com
1-855-845-9302

Specifications subject to change without prior notice.

IMPORTANT ORDERING INSTRUCTIONS

Orders must be placed directly with ROi using ROi reorder numbers to receive HealthTrust contracted pricing. Specifications subject to change without prior notice.



Distributed by
Resource Optimization
& Innovation, LLC
Nashville, TN 37203